



Introduction

- Radiotherapy is a primary, adjuvant or neoadjuvant treatment for a number of different cancers such as glioblastoma, breast, lung and prostate.
- Significant advances have been made in improving the delivery of ionizing radiation to provide precise dosing with reduced side effects to surrounding normal tissue using image guided micro-irradiation (IGMI).
- However in the preclinical setting the use of IGRT is less common with traditional irradiation studies utilising whole body irradiation with lead shielding attempting to focus the radiation to a specific area on the animal or simple single beam techniques.
- The development of the image-guided small animal radiation research platform (SARRP) allows the treatment of animal models of cancer more accurately and importantly, with planned protocols similar to those utilised in the clinic.
- Here we demonstrate the application of IGMI to treat subcutaneous xenograft tumours established from both cell lines and patient-derived material with little or no adverse effects, as well as the utilisation of the SARRP for *in vitro* screening.

SARRP Features

- The SARRP integrates cone beam computed tomography (CBCT) imaging (high resolution, low imaging dose and 3D reconstruction) with radiation treatment (X-ray).
- Irradiation & imaging takes place in a chamber that incorporates a gantry and robotic specimen stage enabling non-coplanar field arrangements and anteriorposterior/posterior-anterior irradiation (Figure 1).
- Image fusion options for easy target localization, dose planning and avoidance of normal organs at risk.
- High precision beam geometry to achieve conformal dose distributions and clinical quality.
- Open platform to enable the addition of other imaging modalities for future research.

Figure 1: External view of SARRP (left) and internal view showing robotic stage, rotating gantry and X-ray tube (right)





Patient relevant preclinical in vivo models using image-guided small animal irradiation for drug discovery

Andrew Mckenzie¹, Nektaria Papadopoulou¹, Lucy Page¹, Jason King¹, Henry Li², Martin Page¹, Ian Wilson³, and Rajendra Kumari^{1*} ¹PRECOS Ltd., a Crown Bioscience Company, Loughborough, UK; ²Crown Bioscience Inc. Santa Clara, US, ³Xstrahl Ltd, Camberley, UK

Studv Da

Methods

•In vitro assay: HCC-827 (NSCLC adenocarcinoma cell line with an activating EGFR mutation, del E746-A750) and HCC827-ER (Erlotinib resistant variant) were grown in T25 flasks and treated with irradiation. Cells were counted after 6 days.


•In vivo xenograft: HCC-827 and HCC-827ER1 cells were implanted subcutaneously in nude mice (HsdOla:MF1-*Foxn1^{nu}*). Tumour measurements and body weights were taken 3 times weekly and treatment initiated when the tumours reached a mean volume of ~200mm³ (n=5 per group). Erlotinib was dosed at 25mg/kg po QD.


•In vivo patient-derived xenograft (PDX): We have established Caucasian non-small cell lung cancer (NSCLC) PDX models which are maintained subcutaneously in vivo in nude mice (HsdOla:MF1-Foxn1^{nu}) admixed with a human stromal cell component (bone marrow-derived human mesenchymal stem cells, ScienCell). Tumour measurements and body weights were taken 3 times weekly and dosing initiated in 2 models when the tumours reached a mean volume of $\sim 200 \text{ mm}^3$.

•In vivo Irradiation: Mice were anaesthetised and transported to the SARRP where CBCT images were acquired. Using the MuriSlice software the isocenter of the tumour was identified and aligned with the central axis of the beam. Fractionated irradiation was administered with the SARRP (225 kV peak X-ray beams; dose rate of 2.5 Gy/min) using collimators of various diameters and a double beam (gantry position at 0° and 180°) under the guidance of the CBCT. A tolerability was performed initially to evaluate 3Gy/ day × 5 days for 2 week.

In Vitro Assay

HCC-827 showed a dose-dependent kill effect in response to elevating levels of irradiation resulting in reduced cell number (Figure 2). Similar the Erlotinib resistant line showed a dose response however the response to 5Gy was markedly less than in the wildt-ype parental line suggesting that the acquired resistance to Erlotinib bestowed some resistance to irradiation.





alone, of in combination with new agents.

*Contact information: Dr. Rajendra Kumari, rajendra.kumari@precos.co.uk

Download the free Layar App