

# 3D Cell Culture Platform

## *In Vitro/Ex Vivo* Assays for Cancer Pharmacology Factsheet



Capture the reality of tumor biology *in situ* by evaluating your compounds in CrownBio's extensive platform of 3D *ex vivo* and *in vitro* assays

Around 95% of new anticancer drugs eventually fail in clinical trials<sup>(1)</sup>, despite robust indications of activity in existing *in vitro* preclinical models, suggesting that innovative models are required that better mimic tumor biology *in vivo*. To increase the speed and efficiency of the preclinical oncology research process a high proportion of drug discovery is currently performed using cell cultures. However, culturing cells on flat plasticware results in artificial 2D monolayering, with cells suffering additional stress as they spread and flatten on the plastic surface of the flask and fail to engage in efficient cell-cell and cell-matrix contacts. All of these factors are known to influence gene expression and protein function therefore raising the question as to whether a 2D cell culture model represents, with sufficient fidelity, the corresponding tumor *in situ*.

Cells that make up tissues follow complex and dynamic 3D arrangements, which are important for their physiology. The 3D architecture of a cell influences its ability to respond to external stimuli and activate the appropriate response. Cancer cells (much like their normal counterparts) benefit from such 3D interactions with each other and with the surrounding extracellular matrix by establishing a unique growing niche, the tumor microenvironment (TME). In order to investigate the efficacy of a new compound as an anticancer agent it is necessary to take in to consideration the 3D nature of a tumor and its communications with the TME.

3D cell cultures represent a significant step forward towards capturing the realities of tumor biology *in situ*. In a 3D cell culture, tumor cells are grown to avoid contact with plastic allowing them to keep their native architecture and to establish interactions with each other, which minimizes stress and an artificial response. A comparison between 2D culturing and the improvement presented by 3D cell culture is shown in **Table 1**.

**Table 1: Comparison of 2D and 3D Cell Culture**

	2D	3D
Cell Shape	Flat. Thickness: 3µm	Ellipsoid. Thickness: 10-30µm
Environment	About 50% of cell surface in contact with culture flask and very small proportion of cell in contact with neighboring cells	Nearly 100% of surface exposed to other cells or to the matrix
Cell Behavior	Differences in differentiation <sup>(2)</sup> , drug metabolism <sup>(3)</sup> , gene expression <sup>(4)</sup> , proliferation rates <sup>(5)</sup> , response to stimuli <sup>(6)</sup> , and viability <sup>(7)</sup>	

### CrownBio 3D Cell Culture Resources

At CrownBio, we offer a variety of *ex vivo* and *in vitro* 3D assay systems for a range of *ex vivo* and *in vitro* models for preclinical oncology research.

Our *ex vivo* models are derived from **HuPrime®** and **HuKemia®** - CrownBio's collection of highly characterized patient-derived xenograft (PDX) models, the largest commercial PDX collection available.

Our models provide a rich resource for the derivation of primary cells and cell lines across a range of cancer types with many molecular targets.

Our *ex vivo* models include freshly isolated cells from PDX tumors (utilized in assays immediately on the day of isolation). We have currently validated almost 70 PDX-derived primary cells for 3D cell culture assays, with cells derived from both solid tumors and blood cancer models (detailed in **Table 2**).

**Table 2: PDX-Derived Primary Cells Validated in 3D Cell Culture Assays**

Models in orange are newly validated.

Tissue of Origin	HuPrime Identifier
ALL	AL5511, AL7443, AL7473
AML	AM5512, AM7440, AM7577
Brain	BN2276, BN2287, BN2331, BN2338
Breast	BR1282, BR1283, BR1458, BR1474
Colorectum	CR0029, CR1520, CR1560, CR1574, CR2161, CR2506, CR2518, CR2520, CR2524, CR2545, CR3085
Esophageal	ES0195, ES0201
Gallbladder	GL0440
Liver	LI0050, LI0574, LI0612, LI0752, LI0941, LI1037, LI1098, LI1646
Lung	LU0038, LU0299, LU0330, LU0367, LU0377, LU0387, LU0395, LU0743, LU0755, LU0884, LU1143, LU1144, LU1147, LU1155, LU1215, LU1271, LU1380, LU1423, LU1868, LU1901, LU2503, LU2505, LU2512, LU3075, LU6409
Ovarian	OV0205, OV0273, OV1658
Pancreas	PA1168
Stomach	GA0037, GA0087, GA0103, GA0114

# 3D Cell Culture Platform *In Vitro/Ex Vivo*

## Assays for Cancer Pharmacology Factsheet

We also provide fresh frozen PDX cells for 3D assay. Cell suspensions are isolated, frozen, and thawed as needed for assays. This can shorten assay timelines to two weeks, with no lag time while PDX tumor donors grow out. CrownBio provides fresh frozen PDX cells from a range of cancer types (shown in **Table 3**), with further models under development.

**Table 3: PDX-Derived Fresh Frozen Cells Validated in 3D Cell Culture Assays**

PDX Model Cancer Type	Model No.	PDX Model Cancer Type	Model No.
Bladder	2	Kidney	1
Brain	2	Liver	11
Breast	5	Lung	20
Cervical	1	Lymphoma	3
Cholangiocarcinoma	3	Metastatic Carcinoma	1
Colorectal	18	Ovarian	6
Esophageal	7	Pancreatic	3
Gallbladder	2	Sarcoma	3
Gastric	10	<b>Total</b>	<b>103</b>
Head and Neck	5		

HuPrime models have also been used to establish our proprietary PrimePanel™ cell line collection<sup>(8)</sup>, consisting of mouse stromal-cell-depleted cancer cell cultures. Our cell lines are all early passage (<10) and maintain essential histopathological features and genetic profiles of the original patient tumors including genomic mutational status, biochemical signaling, and response to tumor cell autonomously targeted therapeutics. PrimePanel provides a versatile, cost-effective, high-throughput screening platform to assess drug efficacy prior to *in vivo* study. At CrownBio we have established five PrimePanel cell lines for 3D use, detailed in **Table 4**.

**Table 4: PrimePanel Cell Lines Validated in 3D Cell Culture Assays**

PrimePanel ID	Cancer Type	Patient and Disease Background
CR5063CL	Colon carcinoma	Female aged 27. Prior treatment with FOLFOX, Camptosar®, Erbitux®, FUDR®, Leucovorin
LU5143CL	NSCLC	Male, aged 78
LU5224CL	SCLC	Female, aged 48. Invasive high grade carcinoma
LU5381CL	NSCLC	Female, aged 57. Prior treatment with IFEX®, Taxol®, Gemzar®
ME5338CL	Melanoma	Metastatic melanoma (brain)

To complement our *ex vivo* models, CrownBio has established 3D cell cultures for >100 commercial cell lines that we have validated for *in vitro* drug efficacy studies, including breast, lung, ovary, and pancreatic cancer cell lines. A summary of our validated *in vitro* 3D cell lines is shown in **Table 5**.

**Table 5: In Vitro Cell Lines Validated in 3D Cell Culture Assays**

	Methylcellulose	Alvetex	Soft Agar
Cancer Type	Cell Line		
Blood	Daudi, DoHH-2, EHEB, Jurkat (clone E6-1), JVM-13, JVM-2, JVM-3, K-562, Karpas299, Kasumi, MEC-1, MEC-2, MEG-01, Molm-16, MOLT-4, Namalwa, OCI-LY19, REH, SU-DHL-10, SU-DHL-5, SUP-B15, TF-1, THP-1, Toledo, U-937, WSU-DLCL-2		
Bone	CADO-ES1, RD-ES		
Brain & Nerves	A172, H4, IMR-32, LN18, LN229, SF126, SH-SY5Y, SK-N-SH, U-118 MG, U251, U-87 MG		
Breast	BT-549, DU4475, MDA-MB-231, MDA-MB-436		MDA-MB-231
Colorectum	HT-29, HCT 116,		HT-29, HCT-116
Esophageal	KYSE70, KYSE270, KYSE410, TE-1		
Head & Neck	FaDu, SW579		
Liver	HCCLM3, Hep G2, HLF, HUH-7, JHH-5, JHH-7, MHCC97H, PLC/PRF/5, SK-HEP-1, SNU-398, SNU-878, SNU-886	Hep 3B	Hep 3B, Hep G2, JHH-7, PLC/PRF/5, SK-HEP-1
Lung	A549, Calu-6, DMS53, DMS79, DMS114, EBC-1, H69AR, HCC4006, NCI-H1155, NCI-H1299, NCI-H1373, NCI-H1395, NCI-H1417, NCI-H1435, NCI-H1437, NCI-H157, NCI-H1581, NCI-H1650, NCI-H1688, NCI-H1703, NCI-H1792, NCI-H1836, NCI-H1975, NCI-H2052, NCI-H209, NCI-2227, NCI-H226, NCI-H23, NCI-H358, NCI-H446, NCI-H460, NCI-H520, NCI-H522, NCI-H526, NCI-H69, NCI-H82, SK-MES-1		NCI-H1299
Ovary	OVCAR-8, SK-OV-3, SW626, SW756		
Pancreas	AsPc-1, Capan-1, MIAPaCa-2, PANC-1, PL45		CFPAC-1, MIAPaCa-2
Prostate	PC-3		
Stomach	A3/KAW, AGS, AZ521, GTL-16, HGC-27, Hs 746T, IM95, IM95m, KATO III, MKN1, MKN45, MKN74, NCI-N87, NUGC-3, NUGC-4, OCUM-1, SCH, SNU-1, SNU-16, SNU-5, SNU-484, SNU601, SNU-620, SNU-638, SNU-668, SNU-719, YCC-2, YCC-7, YCC-10, YCC-11		AGS, AZ-521, GTL-16, HGC-27, MKN1


# 3D Cell Culture Platform *In Vitro/Ex Vivo* Assays for Cancer Pharmacology Factsheet

Our 3D cell culture platforms utilizing the models detailed above for compound screening include the methylcellulose assay, including commercially available systems such as the 3D Alvetex® matrix (which can be used for cells which do not grow well in methylcellulose), soft agar colony formation, and the 3D tumor growth assay (3D TGA).

## Ex Vivo Models in Methylcellulose Assays using HuPrime PDX Models

Our 3D clonogenic *in vitro/ex vivo* methylcellulose assay provides an intimate connection between low cost, high throughput *in vitro* drug screening, and clinically relevant *in vivo* drug efficacy studies with our unique HuPrime models (assay concept shown in **Figure 1**). All of the PDX-derived *ex vivo* models discussed above can be utilized in this assay.





**Figure 1: Ex Vivo Clonogenic Assay Using Freshly Isolated Cells From HuPrime PDX Models**



For example, for freshly isolated primary cells inoculated into immunodeficient mice, xenografts are harvested once the tumor has reached the appropriate size (typically between 400mm<sup>3</sup> and 800mm<sup>3</sup>). The harvested tumor is lysed and cells are seeded on methylcellulose or 3D Alvetex plates at the appropriate dilution to allow growth in 3D, with cell density optimized to provide highly reproducible results. The anticancer agent to be tested can then be administered to cells starting from Day 2 up to Day 10 of culture. Endpoints include phase contrast imaging of cells and CellTiter-Glo® (CTG) assay to verify anticancer compound efficacy.

An example *ex vivo* efficacy study for cisplatin on freshly isolated primary cells from the gallbladder tumor model GL0440 is shown in **Table 6**. Cisplatin was administered to tumor cells at 0.1µM, 1µM, and 100µM and phase contrast images (10X) were taken for each dose to evaluate drug efficacy. Dose response curves can then be generated; **Figure 2** shows the dose response to cisplatin for cells derived from the OV5397 ovarian cancer PDX model in the methylcellulose assay.

**Table 6: Ex Vivo Efficacy Study of Cisplatin on Primary Cells Derived from GL0440 Grown in 3D**

Vehicle Control	Cisplatin 0.1µM	Cisplatin 1µM	Cisplatin 100µM
			


# 3D Cell Culture Platform *In Vitro/Ex Vivo*

## Assays for Cancer Pharmacology Factsheet



**CrownBio**  
CONNECTING SCIENCE TO PATIENTS

**Figure 2: Dose Response Curve for Cisplatin Effect on Primary Cells Derived from OV5397 Grown in 3D**



**Table 7: Comparison of Tumor Conditions Across Humans, 2D, and 3D Culture<sup>(11)</sup>**

Tumor Condition	Human	2D Culture	3D-TME Culture
Oxygen (Hypoxia)	0-5%	21%	≤5%
Stiffness	200-4,000 Pa	3,000,000,000 Pa	200-4,000 Pa
pH	acidic (≤7.0)	buffered (≥7.2)	acidic (≤7.0)
Dimensionality	3D	2D	3D
Glucose availability	<7mM	10-25mM	<7mM


### Ex Vivo Models in 3D TGA Assays using HuPrime PDX models and PrimePanel Cells

3D TGA assays have become a proven platform for oncology drug development, with numerous publications covering research in a range of cancer types<sup>(9,10,11,12)</sup>. The 3D TGA assay utilizes a low stiffness laminin rich extracellular matrix (IrbME, Cultrex®) to embed tumor cells, and admixing with hMSCs (e.g. IL-6, HGF) and CAFs provides the paracrine signaling present in the TME of solid tumors. The addition of hormones (e.g. DHT/E2), restriction of glucose (≤7mM), and maintenance of an acidic pH (6.8) provide a 3D assay that is both “humanized” and TME-aligned for profiling of PDX-derived cells and drug panels<sup>(9)</sup>.

**Table 7** shows a comparison between tumor conditions in humans, 2D culture, and the humanized 3D cell culture environment, with the 3D culture more closely representing the human condition.

Our 3D TGA assay has been optimized to define a linear range for cell growth and for total cell population numbers using alamarBlue® colorimetric assay (further endpoints such as 3D Glo are also available). Individual cell types can also be monitored via fluorescent/bioluminescent read-outs if cells are transduced. The 3D TGA can be utilized with all of the *ex vivo* models described above for drug evaluation and further research as needed. **Figure 3** shows example IC<sub>50</sub> evaluations using our 3D TGA with selected colorectal, lung, and melanoma PrimePanel cell lines.

**Figure 3: Drug Panel IC<sub>50</sub> Evaluation using PrimePanel Cell Lines in a 3D TGA**



PrimePanel Cell Line	IC <sub>50</sub> (μM)										
	5-FU	Carboplatin	Cisplatin	Doxorubicin	Eribitux	Etoposide	Gemcitabine	Irinotecan	Oxaliplatin	Paclitaxel	Vinorelbine
CR5063CL	0.405	-	-	-	~0.157	-	0.020	1.54	~1.99	-	-
LU5224CL	-	~0.014	-	-	-	0.248	0.080	-	-	~0.017	0.021
ME5338CL	4.59	-	>100	8.67	-	~2.09	-	-	-	0.007	-


# 3D Cell Culture Platform *In Vitro/Ex Vivo*

## Assays for Cancer Pharmacology Factsheet

PDX tumor models are ideal for the preclinical assessment of anti-cancer agent efficacy because they closely mimic the pathological phenotypes of the original tumors. However, the labor-intensive, time-consuming, low throughput, and costly nature of *in vivo* animal studies can be limiting for their use in early stage drug development. PrimePanel cells are not generated from all PDX; however, the 3D TGA system enables the use of fresh PDX material in a 96-well format.

Correlation has been shown between response to anticancer agents in 3D TGA and in *in vivo* models and expected clinical outcome. LU6422 is a NSCLC PDX model of the adenocarcinoma subtype. The PDX model harbors a heterozygous c.2573T>G (L858R) in the intracellular kinase domain, which is associated with hypersensitivity to EGFR targeted therapies such as erlotinib in patients. *In vivo*, following subcutaneous implantation with MSCs and treatment with 25 and 50mg/kg erlotinib, complete regression of all tumor growth is observed for this PDX model (**Figure 4A**). Cells derived from LU6422 are also hypersensitive to erlotinib in 3D TGA, correlating with *in vivo* response (**Figure 4B**). These data were recently published with our collaborators in *Molecular Cancer Therapeutics* <sup>(9)</sup>.

**Figure 4: Correlation of LU6422 PDX Response *In Vivo* and in 3D TGA**




### 3D TGA Large Scale Screens

PDX models can be used in mouse clinical trials (HuTrials™) to identify responder and non-responder populations, and to fully elucidate signatures and predictive biomarkers for patient stratification in the clinic. Cells derived from these PDX can also be used in 3D TGA assays, to perform large scale efficacy screens to potentially identify drugs and combinations most likely to benefit individual patients<sup>(9)</sup>.

### *In Vitro* Models in Methylcellulose and Soft Agar Assays using Cell Lines

To complement our *ex vivo* assays, CrownBio has validated a range of commercially available *in vitro* cell lines for 3D assay. Growing cells are harvested and seeded on methylcellulose, soft agar, or other commercially available matrices at the appropriate dilution to achieve a final cell density which is optimized for each cell line (e.g. 5,000 cells per well). Starting from Day 2, anticancer compounds are administered to each well according to desired concentration design. Cell growth is monitored daily until endpoint, when cell number is counted using CTG assay. **Figure 5** shows a dose response curve for the head and neck cancer line FaDu, grown in 3D and treated with increasing doses of an experimental anticancer compound to test its efficacy.

**Figure 5: Dose Response Curve for the FaDu Head and Neck Cancer Cell Line**



We also provide soft agar assays in low-attachment 96 well plates. Cell lines are plated at  $3 \times 10^3$  cells/well in 0.4% agar, over a base layer of 0.6% agar. The test agent is added with the growth media, and endpoints are CTG or imaging assays. **Figure 6** shows comparison images for both endpoints from a soft agar assay, for cisplatin dose response. The calculated EC<sub>50</sub> values are comparable, validating imaging as an endpoint for this assay.




# 3D Cell Culture Platform *In Vitro/Ex Vivo*

## Assays for Cancer Pharmacology Factsheet



**CrownBio**  
CONNECTING SCIENCE TO PATIENTS

**Figure 6: Dose Response Curve for CFPAC-1 Cancer Cell Line Treated with Cisplatin**



### Summary

Tumors are 3D pathological entities. Their 3D nature should be taken into consideration when assessing the efficacy of new anticancer agents as it is able to influence their response to external stimuli such as chemotherapy. Attrition rates in oncology drug development are significantly higher than in other disease areas, with only 5% of anticancer agents which show activity in existing preclinical models going on to be licensed after successful Phase III trials<sup>(1)</sup>. The majority of these failures are related to efficacy rather than toxicity issues, which confirms that current preclinical oncology 2D cell culture systems do not represent, with sufficient fidelity, the corresponding 3D tumor *in situ*.

CrownBio has developed a variety of *ex vivo* and *in vitro* models for 3D culture and preclinical oncology research. Our *ex vivo* models are derived from our well-characterized and validated **HuPrime** and **HuKemia** PDX models, and we provide freshly isolated cell lines, fresh frozen cells, as well as our proprietary **PrimePanel** cell lines which are early passage (<10) and maintain essential histopathological features and genetic profiles of the original patient tumors. To complement these resources, we also provide >100 commercial cell lines which we have validated for 3D culture.

Our *ex vivo* and *in vitro* models can be utilized across a range of 3D assays. CrownBio provides a 3D clonogenic *in vitro/ex vivo* methylcellulose assay which allows an intimate connection between low cost, high throughput *in vitro* drug screening and clinically relevant *in vivo* drug efficacy studies. Soft agar assays provide a similar function, with CellTiter-Glo or imaging endpoints as required.

We also offer the 3D TGA, which has become a proven platform for preclinical oncology drug development. 3D TGA is closely aligned to the TME with epithelial cells being cultured in combination with stromal cells and extra-cellular matrix in serum-free conditions, providing a more clinically relevant culture method. In combination with PDX models, the 3D TGA has the potential to accelerate research by providing a rapid and scalable *ex vivo* drug-sensitivity screen of multiple single and combinatorial agents<sup>(9)</sup>.

CrownBio can be contacted at [busdev@crownbio.com](mailto:busdev@crownbio.com) for any further questions or information on our 3D cell culture platforms or for any other information on other CrownBio products and services.

### References

- <sup>1</sup>Hutchinson L, Kirk R. High drug attrition rates – where are we going wrong? *Nature Reviews Clinical Oncology* 2011;8: 189-190.
- <sup>2</sup>Tian XF, Heng BC, Ge Z *et al.* Comparison of osteogenesis of human embryonic stem cells within 2D and 3D culture systems. *Scandinavian Journal Clinical Laboratory Investigation* 2008;68(1): 58-67.
- <sup>3</sup>Elkayam T, Shaprut S, Dvir-Ginzberg M *et al.* Enhancing drug metabolism activities of C3A – A human hepatocyte cell line - By tissue engineering within scaffolds. *Tissue Engineering* 2006;12(5): 1357-1368
- <sup>4</sup>Shoeters G, Leppens H, Van Gorp U, and Van Den Heuvel R. Haemopoietic long-term bone marrow cultures from adult mice show osteogenic capacity in vitro on 3-dimensional collagen sponges. *Cell Proliferation* 1992;25(6): 587-603.
- <sup>5</sup>Scaglione S, Braccini A, Wendt D *et al.* Engineering of osteoinductive grafts by expansion of ovine bone marrow stromal cells directly on 3D ceramic scaffolds. *Biotechnology and Bioengineering* 2005;93(1): 181-187.
- <sup>6</sup>Debeb BG, Zhang X, Krishnamurthy S *et al.* Characterizing cancer cells with cancer stem cell-like features in 293T human embryonic kidney cells. *Molecular Cancer* 2010;9: 180.
- <sup>7</sup>Laurentin CT, El-Amin SF, Ibim SE *et al.* A highly porous 3-dimensional polyphosphazene polymer matrix for skeletal tissue regeneration. *Journal Biomedical Materials Research* 1996;30(2): 133-138.
- <sup>8</sup>Zhang Y, Ge Y, Liu Y *et al.* PrimePanel provides a high throughput in vitro drug screening platform that intimately links to in vivo pharmacological analysis in PDX models [abstract]. In: Proceedings of the 104th Annual Meeting of the American Association for Cancer Research; 2013 Apr 6-10; Washington, DC. Philadelphia (PA): AACR; *Cancer Research* 2013;73(8 Suppl):Abstract nr 2787.
- <sup>9</sup>Onion D, Argent RH, Reece-Smith AM *et al.* 3-Dimensional Patient-Derived Lung Cancer Assays Reveal Resistance to Standards-of-Care Promoted by Stromal Cells but Sensitivity to Histone Deacetylase Inhibitors. *Molecular Cancer Therapeutics* 2016;15(4): 753-756.
- <sup>10</sup>Casneuf T, Axel AE, King P *et al.* Interleukin-6 is a potential therapeutic target in interleukin-6 dependent, estrogen receptor-α-positive breast cancer. *Breast Cancer: Targets and Therapy* 2016;8: 13–27.
- <sup>11</sup>Krausz E, de Hoogt R, Gustin E *et al.* Translation of a Tumor Microenvironment Mimicking 3D Tumor Growth Co-culture Assay Platform to High-Content Screening. *Journal of Biomolecular Screening* 2013;18(1): 54-66.
- <sup>12</sup>Spaeth EL, Dembinski JL, Sasser AK *et al.* Mesenchymal Stem Cell Transition to Tumor-Associated Fibroblasts Contributes to Fibrovascular Network Expansion and Tumor Progression. *PLoS ONE* 2009;4(4): 34992.

