

In Vitro Potency Assays for Immune Checkpoint Blockade using Human Primary Cells, Murine **HuGEMM™** Immune Cells, and Patient-Derived Tumor Organoids

Abstract
No.237

Xuefei Yan, Hongjuan Zhang, Jun Zhou, Jia Zheng, Shuang Zhu, Yahong Zhang, Rui Zhang, Mingfa Zang, Annie An, Xiaoxi Xu, Kevin Xu, Davy Xuesong Ouyang, Henry Li, and Yujun Huang

Crown Bioscience Inc., 16550 West Bernardo Drive, Building 5, Suite 525, San Diego, CA 92127, USA.

INTRODUCTION

The demand for evaluating the potency of immune checkpoint modulators is steadily growing within immuno-oncology drug development. We aimed to establish a platform to assess the effects of immune checkpoint blockade using human primary immune cells, humanized murine primary immune cells, and co-cultures of tumor cells or patient-derived tumor organoids with immune cells.

METHODS

Antigen-specific T cell recall and activation assay: PBMCs from CMV seropositive and HLA-A02-01+ donors were stimulated with CMVpp65 peptide for one week. CVMpp65-tetramer was used to measure the expansion of CMVpp65-specific CD8⁺ T cells by flow cytometry.

Mixed lymphocyte reaction (MLR) assay: Isolated dendritic cells from donor 1 and T cells from donor 2 were co-cultured to evaluate the activation of allogenic T cells in the presence of anti-PD-1 antibodies. IL-2 and IFN- γ production was measured at different time points by ELISA.


Tumor cell killing assay: PD-L1⁺ tumor cells were co-cultured with allogenic T cells from healthy donors. Anti-PD-1 antibody was used to evaluate the effect of checkpoint blockade by measurement of tumor cell killing. Tumor cells were labeled with CFSE. The percentage of dead and live CFSE⁺ tumor cells were analyzed by flow cytometry.

Dendritic cell (DC) activation assay: Bone marrow from CD40 **HuGEMM** mice were used for differentiation of DCs. DCs were stimulated with an anti-human CD40 agonistic antibody to evaluate the activation and maturation of DC by flow cytometry.

PD-1 blockade assay for co-culture of tumor organoids and PBMCs: PD-L1⁺ luciferase engineered tumor organoids were co-cultured with PBMCs from healthy donors for 48 hours. Luciferase activity was measured to indicate the organoid killing mediated by allogenic T cells in the presence or absence of anti-PD-1. Activation of allogenic cells responding to tumor organoids was also evaluated by IFN- γ ELISA.

RESULTS

Fig 1. PD-1 blockade assay for antigen-specific T cell recall and activation. Anti-PD-1 antibody (Keytruda®) treatment enhances the expansion of CMVpp65-specific CD8⁺ T cells after CMVpp65 peptide stimulation for one week



RESULTS

Fig 2. MLR assay to evaluate PD-1 blockade by anti-PD-1 Abs. (A) IL-2; (B) IFN- γ in the supernatant were measured by ELISA at different time points




Fig 3. Allogenic T-cell mediated tumor cell killing assay to evaluate the effect of PD-1 blockade. CFSE labelled PD-L1⁺ HT-3 tumor cells were co-cultured with PBMC in the presence of anti-PD-1 Ab or isotype for 5 days. % dead cells among CFSE⁺ tumor cells were analyzed with FACS




Fig 4. Bone marrow-derived DC (BMDC) activation assay using **HuGEMM mice.** CD40 **HuGEMM** mice were used to generate BMDC expressing human CD40. Anti-CD40 agonistic antibody (selicrelumab) or LPS was used to stimulate DC and activation markers MHCII, CD80, and CD86 were measured on CD11c⁺ cells





Fig 5. Co-culture assays of tumor organoids and allogenic T cells to evaluate PD-1 blockade. PD-L1⁺ tumor organoids engineered with luciferase were co-cultured with allogenic T cells from different donors. (A) Organoid killing was evaluated by measuring luciferase activity in organoids; (B) Activation of allogenic T cells was assessed by measuring IFN- γ in culture supernatants



SUMMARY

- We have established a panel of *in vitro* assays with human tissue/tumor/targets to evaluate the potency of next generation immune checkpoint inhibitors
- Potency assays for immune checkpoint blockade, such as anti-PD-1 antibodies, were validated using T cell activation and mixed lymphocyte reaction (MLR) assays
- Co-culture of tumor cells or tumor organoids with allogenic T cells was established to measure the effect of PD-1 blockade on tumor cell killing mediated by T cells and activation of T cells
- HuGEMM** mice express engineered human immune checkpoint targets on immune cells and they serve as an excellent resource of primary immune cells to test drug candidates targeting human checkpoints *in vitro*