



TMA006

Mixed CDX Tumor Microarrays for H&E and IHC staining, formalin-fixed and embedded in paraffin.

Tumor Microarray Profile

Catalog #: TMA-CP-MX-006

Model Type: CDX

Host Species: Mouse

Models On TMA: 31

Cores/Model: 3

Cores On TMA: 93

Core Size: 1 mm diameter

Fixation Method: Formalin fixed, paraffin embedded

Applications: H&E and IHC staining

Storage Conditions: 4°C Stable for 12 months from date of purchase when stored as recommended

Slide Preparation: Slides are sealed with wax; heat slide at 60° for 1 hour before de-waxing

QC: Each TMA slide has >90% core occupancy. H&E and Ki67 IHC staining have been performed on a representative slide of each lot of TMA product.

Tumor Microarray Map

Slide Label													
	TMA-CP-MX-006, Human				triplicate								
	A	B	C	D	E	F	G	H	I	J	K	L	
1	KYSE-270 HN/Esophagus SubQ	KYSE410 HN/Esophagus SubQ	OE19 HN/Esophagus/ Gastric cardia SubQ	C666-1 HN/Nose SubQ	Colo-205 Large intestine/Colore ctum SubQ	HCT-15 Large intestine/Colore ctum SubQ	HT-29 Large intestine/Colore ctum SubQ	LS411N Large intestine/Colore ctum SubQ	SNU-81 Large intestine/Colore ctum SubQ	SW480 Large intestine/Colore ctum SubQ	SW948 Large intestine/Colore ctum SubQ	KM12/L4 Colon Cancer SubQ	
	Lovo Colon Cancer SubQ	Hep3B Liver Cancer SubQ	HepG2/C3A Liver Cancer SubQ	JHH-5 Liver Cancer SubQ	MHCC97H Liver Cancer SubQ	PLC/PRF/5 Liver Cancer SubQ	SNU-878 Liver Cancer SubQ	HuCCT1 Liver/Bile duct Cancer SubQ	HPAF-II Pancreas Cancer SubQ	KP4 Pancreas Cancer SubQ	Panc 04.03 Pancreas Cancer SubQ	PANC-1 Pancreas Cancer Orth	
	Capan-2 Pancreas Cancer SubQ	BX-PC-3 Pancreas Cancer SubQ	Capan-1 Pancreas Cancer SubQ	NCI-N87 Stomach Cancer SubQ	SNU-16 Stomach Cancer SubQ	SNU-5 Stomach Cancer SubQ	SNU-620 Stomach Cancer SubQ	KYSE-270 HN/Esophagu s SubQ	KYSE410 HN/Esophagu s SubQ	OE19 HN/Esophagu s/Gastric cardia SubQ	C666-1 HN/Nose SubQ	Colo-205 Large intestine/Color ectum SubQ	
	HCT-15 Large intestine/Color ectum SubQ	HT-29 Large intestine/Color ectum SubQ	LS411N Large intestine/Color ectum SubQ	SNU-81 Large intestine/Color ectum SubQ	SW480 Large intestine/Color ectum SubQ	SW948 Large intestine/Color ectum SubQ	KM12/L4 Colon Cancer SubQ	Lovo Colon Cancer SubQ	Hep3B Liver Cancer SubQ	HepG2/C3A Liver Cancer SubQ	JHH-5 Liver Cancer SubQ	MHCC97H Liver Cancer SubQ	
	PLC/PRF/5 Liver Cancer SubQ	SNU-878 Liver Cancer SubQ	HuCCT1 Liver/Bile duct Cancer SubQ	HPAF-II Pancreas Cancer SubQ	KP4 Pancreas Cancer SubQ	Panc 04.03 Pancreas Cancer SubQ	PANC-1 Pancreas Cancer Orth	Capan-2 Pancreas Cancer SubQ	BX-PC-3 Pancreas Cancer SubQ	Capan-1 Pancreas Cancer SubQ	NCI-N87 Stomach Cancer SubQ	SNU-16 Stomach Cancer SubQ	
	SNU-5 Stomach Cancer SubQ	SNU-620 Stomach Cancer SubQ	KYSE-270 HN/Esophagus SubQ	KYSE410 HN/Esophagus SubQ	OE19 HN/Esophagus/ Gastric cardia SubQ	C666-1 HN/Nose SubQ	Colo-205 Large intestine/Colore ctum SubQ	HCT-15 Large intestine/Colore ctum SubQ	HT-29 Large intestine/Colore ctum SubQ	LS411N Large intestine/Colore ctum SubQ	SNU-81 Large intestine/Colore ctum SubQ	SW480 Large intestine/Colore ctum SubQ	
	SW948 Large intestine/Colore ctum SubQ	KM12/L4 Colon Cancer SubQ	Lovo Colon Cancer SubQ	Hep3B Liver Cancer SubQ	HepG2/C3A Liver Cancer SubQ	JHH-5 Liver Cancer SubQ	MHCC97H Liver Cancer SubQ	PLC/PRF/5 Liver Cancer SubQ	SNU-878 Liver Cancer SubQ	HuCCT1 Liver/Bile duct Cancer SubQ	HPAF-II Pancreas Cancer SubQ	KP4 Pancreas Cancer SubQ	
	Panc 04.03 Pancreas Cancer SubQ	PANC-1 Pancreas Cancer Orth	Capan-2 Pancreas Cancer SubQ	BX-PC-3 Pancreas Cancer SubQ	Capan-1 Pancreas Cancer SubQ	NCI-N87 Stomach Cancer SubQ	SNU-16 Stomach Cancer SubQ	SNU-5 Stomach Cancer SubQ	SNU-620 Stomach Cancer SubQ				
8													Mark point

Map of mixed cancer types tumor microarray slide

Please follow standard IHC staining procedure.

Visit <https://www.crownbio.com/databases/xenobase> to view detailed information for each model. Access to XenoBase® requires online registration via the link provided.