



TMA008

Mixed CDX Tumor Microarrays for H&E and IHC staining, formalin-fixed and embedded in paraffin.

Tumor Microarray Profile

Catalog #: TMA-CP-MX-008

Model Type: CDX

Host Species: Mouse

Models On TMA: 28

Cores/Model: 3

Cores On TMA: 87

Core Size: 1 mm diameter

Fixation Method: Formalin fixed, paraffin embedded

Applications: H&E and IHC staining

Storage Conditions: 4°C Stable for 12 months from date of purchase when stored as recommended

Slide Preparation: Slides are sealed with wax; heat slide at 60° for 1 hour before de-waxing

QC: Each TMA slide has >90% core occupancy. H&E and Ki67 IHC staining have been performed on a representative slide of each lot of TMA product.

Tumor Microarray Map

Slide Label													
	TMA-CP-MX-008, Human			triplicate									
	A	B	C	D	E	F	G	H	I	J	K	L	
1	CA46 Blood/Lymphoma SubQ	Caki-1 Burkitt lymphoma SubQ	DOHH-2 Lymphoma SubQ	Granta-519 Blood/Lymphoma SubQ	HL-60 Leukemia SubQ	HuT 102 Blood/Lymphoma SubQ	Karpas299 Lymphoma SubQ	KARPAS-422 Blood/Lymphoma SubQ	Kasumi-1 Leukemia SubQ	MM1.R Blood/Myeloma SubQ	MOLM-13 Blood/Leukemia SubQ	MOLM-16 Leukemia SubQ	
2	MOLP8 Blood/Myeloma SubQ	Molt-4 Blood/Leukemia SubQ	MV-4-11 Leukemia SubQ	NCI-H929 Blood/Myeloma SubQ	OCI-AML3 Blood/Leukemia SubQ	OCI-LY19 Blood/Lymphoma SubQ	OCI-LY7 Blood/Lymphoma SubQ	NA	Pfeiffer Blood/Lymphoma SubQ	Raji Lymphoma SubQ	RPMI-8226 Myeloma SubQ	ST486 Blood/Lymphoma SubQ	
3	SU-DHL-10 Blood/Lymphoma SubQ	THP-1 Blood/Leukemia SubQ	U-937 Lymphoma SubQ	WSU-DLCL2 Blood/Lymphoma SubQ	JeKo-1 Lymphoma SubQ	CA46 Blood/Lymphoma SubQ	Caki-1 Burkitt lymphoma SubQ	DOHH-2 Lymphoma SubQ	Granta-519 Blood/Lymphoma SubQ	HL-60 Leukemia SubQ	HuT 102 Blood/Lymphoma SubQ	Karpas299 Lymphoma SubQ	
4	KARPAS-422 Blood/Lymphoma SubQ	Kasumi-1 Leukemia SubQ	MM1.R Blood/Myeloma SubQ	MOLM-13 Blood/Leukemia SubQ	MOLM-16 Leukemia SubQ	MOLP8 Blood/Myeloma SubQ	Molt-4 Blood/Leukemia SubQ	MV-4-11 Leukemia SubQ	NCI-H929 Blood/Myeloma SubQ	OCI-AML3 Blood/Leukemia SubQ	OCI-LY19 Blood/Lymphoma SubQ	OCI-LY7 Blood/Lymphoma SubQ	
5	NA	Pfeiffer Blood/Lymphoma SubQ	Raji Lymphoma SubQ	RPMI-8226 Myeloma SubQ	ST486 Blood/Lymphoma SubQ	SU-DHL-10 Blood/Lymphoma SubQ	THP-1 Blood/Leukemia SubQ	U-937 Lymphoma SubQ	WSU-DLCL2 Blood/Lymphoma SubQ	JeKo-1 Lymphoma SubQ	CA46 Blood/Lymphoma SubQ	Caki-1 Burkitt lymphoma SubQ	
6	DOHH-2 Lymphoma SubQ	Granta-519 Blood/Lymphoma SubQ	HL-60 Leukemia SubQ	HuT 102 Blood/Lymphoma SubQ	Karpas299 Lymphoma SubQ	KARPAS-422 Blood/Lymphoma SubQ	Kasumi-1 Leukemia SubQ	MM1.R Blood/Myeloma SubQ	MOLM-13 Blood/Leukemia SubQ	MOLM-16 Leukemia SubQ	MOLP8 Blood/Myeloma SubQ	Molt-4 Blood/Leukemia SubQ	
7	MV-4-11 Leukemia SubQ	NCI-H929 Blood/Myeloma SubQ	OCI-AML3 Blood/Leukemia SubQ	OCI-LY19 Blood/Lymphoma SubQ	OCI-LY7 Blood/Lymphoma SubQ	NA	Pfeiffer Blood/Lymphoma SubQ	Raji Lymphoma SubQ	RPMI-8226 Myeloma SubQ	ST486 Blood/Lymphoma SubQ	SU-DHL-10 Blood/Lymphoma SubQ	THP-1 Blood/Leukemia SubQ	
8	U-937 Lymphoma SubQ	WSU-DLCL2 Blood/Lymphoma SubQ	JeKo-1 Lymphoma SubQ										
													Mark point

Map of mixed cancer types tumor microarray slide

Please follow standard IHC staining procedure.

Visit <https://www.crownbio.com/databases/xenobase> to view detailed information for each model. Access to XenoBase® requires online registration via the link provided.