



TMA-CP-MX-009

Mixed CDX Tumor Microarrays for H&E and IHC staining, formalin-fixed and embedded in paraffin.

Tumor Microarray Profile

Catalog #: TMA-CP-MX-009

Model Type: CDX

Host Species: Mouse

Models On TMA: 32

Cores/Model: 3

Cores On TMA: 96

Core Size: 1 mm diameter

Fixation Method: Formalin fixed, paraffin embedded

Applications: H&E and IHC staining

Storage Conditions: 4°C Stable for 12 months from date of purchase when stored as recommended

Slide Preparation: Slides are sealed with wax; heat slide at 60° for 1 hour before de-waxing

QC: Each TMA slide has >90% core occupancy. H&E and Ki67 IHC staining have been performed on a representative slide of each lot of TMA product.

Tumor Microarray Map

Slide Label													
		TMA-CP-MX-009, Human			triplicate								
	A	B	C	D	E	F	G	H	I	J	K	L	
1	Colo-205 Intestine/Large/Colorectum SubQ	HCT-15 Intestine/Large/Colorectum SubQ	HT-29 Intestine/Large/Colorectum SubQ	SNU-81 Intestine/Large/Colorectum SubQ	SW948 Intestine/Large/Colorectum SubQ	Lovo Intestine/Large/Colorectum SubQ	LS411N Intestine/Large/Colorectum SubQ	COLO 320DM Intestine/Large/Colorectum SubQ	GP2d Intestine/Large/Colorectum SubQ	LS 513 Intestine/Large/Colorectum SubQ	SW-48 Intestine/Large/Colorectum SubQ	KYSE-270 Esophagus SubQ	
2	KYSE410 Esophagus SubQ	OE19 Esophagus SubQ	NCI-N87 Gastric SubQ	SNU-5 Gastric SubQ	SNU-620 Gastric SubQ	FaDu HN/Pharynx SubQ	8305 C HN/Thyroid SubQ	Hep3B Liver SubQ	HepG2/C3A Liver SubQ	JHH-5 Liver SubQ	PLC/PRF/5 Liver SubQ	SNU-878 Liver SubQ	
3	HuCCT1 Liver/Cholangiocarcinoma SubQ	HPAF-II Pancreas SubQ	KP4 Pancreas SubQ	Panc 04.03 Pancreas Orth	PANC-1 Pancreas SubQ	BX-PC-3 Pancreas SubQ	Capan-1 Pancreas SubQ	Capan-2 Pancreas SubQ	Colo-205 Intestine/Large/Colorectum SubQ	HCT-15 Intestine/Large/Colorectum SubQ	HT-29 Intestine/Large/Colorectum SubQ	SNU-81 Intestine/Large/Colorectum SubQ	
4	SW948 Intestine/Large/Colorectum SubQ	Lovo Intestine/Large/Colorectum SubQ	LS411N Intestine/Large/Colorectum SubQ	COLO 320DM Intestine/Large/Colorectum SubQ	GP2d Intestine/Large/Colorectum SubQ	LS 513 Intestine/Large/Colorectum SubQ	SW-48 Intestine/Large/Colorectum SubQ	KYSE-270 Esophagus SubQ	KYSE410 Esophagus SubQ	OE19 Esophagus SubQ	NCI-N87 Gastric SubQ	SNU-5 Gastric SubQ	
5	SNU-620 Gastric SubQ	FaDu HN/Pharynx SubQ	8305 C HN/Thyroid SubQ	Hep3B Liver SubQ	HepG2/C3A Liver SubQ	JHH-5 Liver SubQ	PLC/PRF/5 Liver SubQ	SNU-878 Liver SubQ	HuCCT1 Liver/Cholangiocarcinoma SubQ	HPAF-II Pancreas SubQ	KP4 Pancreas SubQ	Panc 04.03 Pancreas Orth	
6	PANC-1 Pancreas SubQ	BX-PC-3 Pancreas SubQ	Capan-1 Pancreas SubQ	Capan-2 Pancreas SubQ	Colo-205 Intestine/Large/Colorectum SubQ	HCT-15 Intestine/Large/Colorectum SubQ	HT-29 Intestine/Large/Colorectum SubQ	SNU-81 Intestine/Large/Colorectum SubQ	SW948 Intestine/Large/Colorectum SubQ	Lovo Intestine/Large/Colorectum SubQ	LS411N Intestine/Large/Colorectum SubQ	COLO 320DM Intestine/Large/Colorectum SubQ	
7	GP2d Intestine/Large/Colorectum SubQ	LS 513 Intestine/Large/Colorectum SubQ	SW-48 Intestine/Large/Colorectum SubQ	KYSE-270 Esophagus SubQ	KYSE410 Esophagus SubQ	OE19 Esophagus SubQ	NCI-N87 Gastric SubQ	SNU-5 Gastric SubQ	SNU-620 Gastric SubQ	FaDu HN/Pharynx SubQ	8305 C HN/Thyroid SubQ	Hep3B Liver SubQ	
8	HepG2/C3A Liver SubQ	JHH-5 Liver SubQ	PLC/PRF/5 Liver SubQ	SNU-878 Liver SubQ	HuCCT1 Liver/Cholangiocarcinoma SubQ	HPAF-II Pancreas SubQ	KP4 Pancreas SubQ	Panc 04.03 Pancreas Orth	PANC-1 Pancreas SubQ	BX-PC-3 Pancreas SubQ	Capan-1 Pancreas SubQ	Capan-2 Pancreas SubQ	
													Mark point

Map of mixed cancer types tumor microarray slide

Please follow standard IHC staining procedure.

Visit <https://www.crownbio.com/databases/xenobase> to view detailed information for each model. Access to

XenoBase® requires online registration via the link provided.