



TMA046

Bladder Cancer Tumor Microarrays derived from PDX for H&E and IHC staining, formalin-fixed and embedded in paraffin.

Tumor Microarray Profile

Catalog #: TMA-HP-BL-046

Model Type: PDX

Host Species: Mouse

Models On TMA: 24

Cores/Model: 3

Cores On TMA: 72

Core Size: 1 mm diameter

Fixation Method: Formalin fixed, paraffin embedded

Applications: H&E and IHC staining

Storage Conditions: 4°C Stable for 12 months from date of purchase when stored as recommended

Slide Preparation: Slides are sealed with wax; heat slide at 60° for 1 hour before de-waxing

QC: Each TMA slide has >90% core occupancy. H&E and Ki67 IHC staining have been performed on a representative slide of each lot of TMA product.

Tumor Microarray Map

Slide Label		TMA-HP-BL-046	triplicate										
		A	B	C	D	E	F	G	H	I	J	K	L
1	BL0597-P8 Bladder cancer	BL2350-P5 Bladder cancer	BL2393-P4 Bladder cancer	BL3325-P7 Bladder cancer	BL3338-P11 Bladder cancer	BL3479-P5 Bladder cancer	BL3578-P4 Bladder cancer	BL5001-P6 Bladder cancer	BL5002-P7 Bladder cancer	BL5003-P10 Bladder cancer	BL5007-P6 Bladder cancer	BL5384-P7 Bladder cancer	
2	BL5422-P6 Bladder cancer	BL5428-P6 Bladder cancer	BL9200-P5 Bladder cancer	BL9201-P3 Bladder cancer	BL9202-P2 Bladder cancer	BL9206-P2 Bladder cancer	BL9208-P4 Bladder cancer	BL9211-P5 Bladder cancer	BL9213-P5 Bladder cancer	BL9214-P3 Bladder cancer	BL9215-P6 Bladder cancer	BL9217-P4 Bladder cancer	
3	BL0597-P8 Bladder cancer	BL2350-P5 Bladder cancer	BL2393-P4 Bladder cancer	BL3325-P7 Bladder cancer	BL3338-P11 Bladder cancer	BL3479-P5 Bladder cancer	BL3578-P4 Bladder cancer	BL5001-P6 Bladder cancer	BL5002-P7 Bladder cancer	BL5003-P10 Bladder cancer	BL5007-P6 Bladder cancer	BL5384-P7 Bladder cancer	
4	BL5422-P6 Bladder cancer	BL5428-P6 Bladder cancer	BL9200-P5 Bladder cancer	BL9201-P3 Bladder cancer	BL9202-P2 Bladder cancer	BL9206-P2 Bladder cancer	BL9208-P4 Bladder cancer	BL9211-P5 Bladder cancer	BL9213-P5 Bladder cancer	BL9214-P3 Bladder cancer	BL9215-P6 Bladder cancer	BL9217-P4 Bladder cancer	
5	BL0597-P8 Bladder cancer	BL2350-P5 Bladder cancer	BL2393-P4 Bladder cancer	BL3325-P7 Bladder cancer	BL3338-P11 Bladder cancer	BL3479-P5 Bladder cancer	BL3578-P4 Bladder cancer	BL5001-P6 Bladder cancer	BL5002-P7 Bladder cancer	BL5003-P10 Bladder cancer	BL5007-P6 Bladder cancer	BL5384-P7 Bladder cancer	
6	BL5422-P6 Bladder cancer	BL5428-P6 Bladder cancer	BL9200-P5 Bladder cancer	BL9201-P3 Bladder cancer	BL9202-P2 Bladder cancer	BL9206-P2 Bladder cancer	BL9208-P4 Bladder cancer	BL9211-P5 Bladder cancer	BL9213-P5 Bladder cancer	BL9214-P3 Bladder cancer	BL9215-P6 Bladder cancer	BL9217-P4 Bladder cancer	
													Mark point

Map of bladder cancer tumor microarray

Please follow standard IHC staining procedure.

Visit <https://www.crownbio.com/databases/hubase> to view detailed information for each model. Access to HuBase™ requires online registration via the link provided.