



TMA051

Brain Cancer Microarrays derived from PDX for H&E and IHC staining, formalin-fixed and embedded in paraffin.

Tumor Microarray Profile

Catalog #: TMA-HP-BN-051

Model Type: PDX

Host Species: Mouse

Models On TMA: 30

Cores/Model: 3

Cores On TMA: 96

Core Size: 1 mm diameter

Fixation Method: Formalin fixed, paraffin embedded

Applications: H&E and IHC staining

Storage Conditions: 4°C Stable for 12 months from date of purchase when stored as recommended

Slide Preparation: Slides are sealed with wax; heat slide at 60° for 1 hour before de-waxing

QC: Each TMA slide has >90% core occupancy. H&E and Ki67 IHC staining have been performed on a representative slide of each lot of TMA product.

Tumor Microarray Map

Slide Label		TMA-HP-BN-051		triplicate									
		A	B	C	D	E	F	G	H	I	J	K	L
	1	BN0769-P9 Brain cancer	BN2276-P11 Brain cancer	BN2287-P16 Brain cancer	BN2289-P27 Brain cancer	BN2331-P13 Brain cancer	BN2338-P11 Brain cancer	BN3733-P2 Brain cancer	BN3801-P8 Brain cancer	BN4027-P8 Brain cancer	BN4030-P10 Brain cancer	BN4043-P3 Brain cancer	BN5463-P5 Brain cancer
	2	NA	BN9222-P3 Brain cancer	BN9224-P4 Brain cancer	BN10369-P2 Brain cancer	BN10394-P3 Brain cancer	BN10401-P3 Brain cancer	BN10410-P3 Brain cancer	BN10431-P2 Brain cancer	BN10438-P2 Brain cancer	BN10475-P3 Brain cancer	BN13331-P2 Brain cancer	BN14976-P5 Brain cancer
	3	BN0769-P9 Brain cancer	BN2276-P11 Brain cancer	BN2287-P16 Brain cancer	BN2289-P27 Brain cancer	BN2331-P13 Brain cancer	BN2338-P11 Brain cancer	BN3733-P2 Brain cancer	BN3801-P8 Brain cancer	BN4027-P8 Brain cancer	BN4030-P10 Brain cancer	BN4043-P3 Brain cancer	BN5463-P5 Brain cancer
	4	NA	BN9222-P3 Brain cancer	BN9224-P4 Brain cancer	BN10369-P2 Brain cancer	BN10394-P3 Brain cancer	BN10401-P3 Brain cancer	BN10410-P3 Brain cancer	BN10431-P2 Brain cancer	BN10438-P2 Brain cancer	BN10475-P3 Brain cancer	BN13331-P2 Brain cancer	BN14976-P5 Brain cancer
	5	BN0769-P9 Brain cancer	BN2276-P11 Brain cancer	BN2287-P16 Brain cancer	BN2289-P27 Brain cancer	BN2331-P13 Brain cancer	BN2338-P11 Brain cancer	BN3733-P2 Brain cancer	BN3801-P8 Brain cancer	BN4027-P8 Brain cancer	BN4030-P10 Brain cancer	BN4043-P3 Brain cancer	BN5463-P5 Brain cancer
	6	NA	BN9222-P3 Brain cancer	BN9224-P4 Brain cancer	BN10369-P2 Brain cancer	BN10394-P3 Brain cancer	BN10401-P3 Brain cancer	BN10410-P3 Brain cancer	BN10431-P2 Brain cancer	BN10438-P2 Brain cancer	BN10475-P3 Brain cancer	BN13331-P2 Brain cancer	BN14976-P5 Brain cancer
													Mark point

Map of brain cancer tumor microarray

Please follow standard IHC staining procedure.

Visit <https://www.crownbio.com/databases> to view detailed information for each model. Access to the databases requires online registration via the link provided.