



TMA069

Breast Cancer Tumor Microarrays derived from PDX for H&E and IHC staining, formalin-fixed and embedded in paraffin.

Tumor Microarray Profile

Catalog #: TMA-HP-BR-069

Model Type: PDX

Host Species: Mouse

Models On TMA: 32

Cores/Model: 3

Cores On TMA: 60

Core Size: 1.5 mm diameter

Fixation Method: Formalin fixed, paraffin embedded

Applications: H&E and IHC staining

Storage Conditions: 4°C Stable for 12 months from date of purchase when stored as recommended

Slide Preparation: Slides are sealed with wax; heat slide at 60° for 1 hour before de-waxing

QC: Each TMA slide has >90% core occupancy. H&E and Ki67 IHC staining have been performed on a representative slide of each lot of TMA product.

Tumor Microarray Map

Slide Label	TMA-HP-BR-069			triplicate									
	A	B	C	D	E	F	G	H	I	J	K	L	
1	BR10007-P5 Breast cancer	BR10553-P5 Breast cancer	BR10596-P3 Breast cancer	BR10670-P5 Breast cancer	BR10679-P6 Breast cancer	BR10688-P5 Breast cancer	BR10698-P5 Breast cancer	BR10699-P4 Breast cancer	BR10703-P5 Breast cancer	BR10749-P4 Breast cancer	BR10767-P5 Breast cancer	BR5009-P7 Breast cancer	
2	BR5010-P13 Breast cancer	BR5011-P10 Breast cancer	BR5017-P13 Breast cancer	BR5020-P8 Breast cancer	BR5023-P5 Breast cancer	BR5337-P6 Breast cancer	BR5380-P10 Breast cancer	BR5400-P7 Breast cancer	BR5494-P13 Breast cancer	BR9225-P7 Breast cancer	BR9227-P6 Breast cancer	BR9236-P8 Breast cancer	
3	BR9462-P10 Breast cancer	BR9473-P5 Breast cancer	BR9493-P20 Breast cancer	BR9552-P4 Breast cancer	BR9565-P4 Breast cancer	BR9655-P3 Breast cancer	BR9690-P9 Breast cancer	BR9801-P6 Breast cancer	BR10007-P5 Breast cancer	BR10553-P5 Breast cancer	BR10596-P3 Breast cancer	BR10670-P5 Breast cancer	
4	BR10679-P6 Breast cancer	BR10688-P5 Breast cancer	BR10698-P5 Breast cancer	BR10699-P4 Breast cancer	BR10703-P5 Breast cancer	BR10749-P4 Breast cancer	BR10767-P5 Breast cancer	BR5009-P7 Breast cancer	BR5010-P13 Breast cancer	BR5011-P10 Breast cancer	BR5017-P13 Breast cancer	BR5020-P8 Breast cancer	
5	BR5023-P5 Breast cancer	BR5337-P6 Breast cancer	BR5380-P10 Breast cancer	BR5400-P7 Breast cancer	BR5494-P13 Breast cancer	BR9225-P7 Breast cancer	BR9227-P6 Breast cancer	BR9236-P8 Breast cancer	BR9462-P10 Breast cancer	BR9473-P5 Breast cancer	BR9493-P20 Breast cancer	BR9552-P4 Breast cancer	
6	BR9565-P4 Breast cancer	BR9655-P3 Breast cancer	BR9690-P9 Breast cancer	BR9801-P6 Breast cancer	BR10007-P5 Breast cancer	BR10553-P5 Breast cancer	BR10596-P3 Breast cancer	BR10670-P5 Breast cancer	BR10679-P6 Breast cancer	BR10688-P5 HP-F4047-2	BR10698-P5 Breast cancer	BR10699-P4 Breast cancer	
7	BR10703-P5 Breast cancer	BR10749-P4 Breast cancer	BR10767-P5 Breast cancer	BR5009-P7 Breast cancer	BR5010-P13 Breast cancer	BR5011-P10 Breast cancer	BR5017-P13 Breast cancer	BR5020-P8 Breast cancer	BR5023-P5 Breast cancer	BR5337-P6 Breast cancer	BR5380-P10 Breast cancer	BR5400-P7 Breast cancer	
8	BR5494-P13 Breast cancer	BR9225-P7 Breast cancer	BR9227-P6 Breast cancer	BR9236-P8 Breast cancer	BR9462-P10 Breast cancer	BR9473-P5 Breast cancer	BR9493-P20 Breast cancer	BR9552-P4 Breast cancer	BR9565-P4 Breast cancer	BR9655-P3 Breast cancer	BR9690-P9 Breast cancer	BR9801-P6 Breast cancer	
													Mark point

Map of breast cancer tumor microarray

Please follow standard IHC staining procedure.

Visit <https://www.crownbio.com/databases/hubase> to view detailed information for each model. Access to HuBase™ requires online registration via the link provided.