



TMA065

Esophageal Cancer Microarrays derived from PDX for H&E and IHC staining, formalin-fixed and embedded in paraffin.

Tumor Microarray Profile

Catalog #: TMA-HP-ES-065

Model Type: PDX

Host Species: Mouse

Models On TMA: 32

Cores/Model: 3

Cores On TMA: 96

Core Size: 1 mm diameter

Fixation Method: Formalin fixed, paraffin embedded

Applications: H&E and IHC staining

Storage Conditions: 4°C Stable for 12 months from date of purchase when stored as recommended

Slide Preparation: Slides are sealed with wax; heat slide at 60° for 1 hour before de-waxing

QC: Each TMA slide has >90% core occupancy. H&E and Ki67 IHC staining have been performed on a representative slide of each lot of TMA product.

Tumor Microarray Map

Slide Label													
		TMA-HP-ES-065		triplicate									
		A	B	C	D	E	F	G	H	I	J	K	L
1	ES0026-P9 Esophageal cancer	ES0042-P17 Esophageal cancer	ES0110-P7 Esophageal cancer	ES0136-P6 Esophageal cancer	ES0148-P10 Esophageal cancer	ES0168-P8 Esophageal cancer	ES0176-P10 Esophageal cancer	ES0178-P11 Esophageal cancer	ES0190-P6 Esophageal cancer	ES0195-P9 Esophageal cancer	ES0199-P6 Esophageal cancer	ES0201-P8 Esophageal cancer	
2	ES0204-P18 Esophageal cancer	ES0214-P11 Esophageal cancer	ES0218-P10 Esophageal cancer	ES0219-P10 Esophageal cancer	ES0630-P9 Esophageal cancer	ES2263-P11 Esophageal cancer	ES2267-P11 Esophageal cancer	ES2311-P10 Esophageal cancer	ES2356-P8 Esophageal cancer	ES3431-P9 Esophageal cancer	ES3862-P11 Esophageal cancer	ES3865-P13 Esophageal cancer	
3	ES3870-P3 Esophageal cancer	ES3897-P8 Esophageal cancer	ES3904-P7 Esophageal cancer	ES3905-P6 Esophageal cancer	ES3910-P8 Esophageal cancer	ES6463-P25 Esophageal cancer	ES6824-P11 Esophageal cancer	ES11086-P2 Esophageal cancer	ES0026-P9 Esophageal cancer	ES0042-P17 Esophageal cancer	ES0110-P7 Esophageal cancer	ES0136-P6 Esophageal cancer	
4	ES0148-P10 Esophageal cancer	ES0168-P8 Esophageal cancer	ES0176-P10 Esophageal cancer	ES0178-P11 Esophageal cancer	ES0190-P6 Esophageal cancer	ES0195-P9 Esophageal cancer	ES0199-P6 Esophageal cancer	ES0201-P8 Esophageal cancer	ES0204-P18 Esophageal cancer	ES0214-P11 Esophageal cancer	ES0218-P10 Esophageal cancer	ES0219-P10 Esophageal cancer	
5	ES0630-P9 Esophageal cancer	ES2263-P11 Esophageal cancer	ES2267-P11 Esophageal cancer	ES2311-P10 Esophageal cancer	ES2356-P8 Esophageal cancer	ES3431-P9 Esophageal cancer	ES3862-P11 Esophageal cancer	ES3865-P13 Esophageal cancer	ES3870-P3 Esophageal cancer	ES3897-P8 Esophageal cancer	ES3904-P7 Esophageal cancer	ES3905-P6 Esophageal cancer	
6	ES3910-P8 Esophageal cancer	ES6463-P25 Esophageal cancer	ES6824-P11 Esophageal cancer	ES11086-P2 Esophageal cancer	ES0026-P9 Esophageal cancer	ES0042-P17 Esophageal cancer	ES0110-P7 Esophageal cancer	ES0136-P6 Esophageal cancer	ES0148-P10 Esophageal cancer	ES0168-P8 Esophageal cancer	ES0176-P10 Esophageal cancer	ES0178-P11 Esophageal cancer	
7	ES0190-P6 Esophageal cancer	ES0195-P9 Esophageal cancer	ES0199-P6 Esophageal cancer	ES0201-P8 Esophageal cancer	ES0204-P18 Esophageal cancer	ES0214-P11 Esophageal cancer	ES0218-P10 Esophageal cancer	ES0219-P10 Esophageal cancer	ES0630-P9 Esophageal cancer	ES2263-P11 Esophageal cancer	ES2267-P11 Esophageal cancer	ES2311-P10 Esophageal cancer	
8	ES2356-P8 Esophageal cancer	ES3431-P9 Esophageal cancer	ES3862-P11 Esophageal cancer	ES3865-P13 Esophageal cancer	ES3870-P3 Esophageal cancer	ES3897-P8 Esophageal cancer	ES3904-P7 Esophageal cancer	ES3905-P6 Esophageal cancer	ES3910-P8 Esophageal cancer	ES6463-P25 Esophageal cancer	ES6824-P11 Esophageal cancer	ES11086-P2 Esophageal cancer	
													Mark point

Map of esophageal cancer tumor microarray

Please follow standard IHC staining procedure.

Visit <https://www.crownbio.com/databases> to view detailed information for each model. Access to the databases requires online registration via the link provided.