



TMA-HP-LI-063

This liver cancer tumor microarray is derived from patient derived xenograft (PDX) tumor samples. Each tumor sample is formalin fixed, paraffin embedded, and spotted in triplicate on the TMA.

Tumor Microarray Profile

Catalog #: TMA-HP-LI-063

Model Type: PDX

Host Species: Mouse

Models On TMA: 32

Cores/Model: 3

Cores On TMA: 96

Core Size: 1 mm diameter

Fixation Method: Formalin fixed, paraffin embedded

Applications: H&E and IHC staining

Storage Conditions: 4°C Stable for 12 months from date of purchase when stored as recommended

Slide Preparation: Slides are sealed with wax; heat slide at 60° for 1 hour before de-waxing

QC: Each TMA slide has >90% core occupancy. H&E and Ki67 IHC staining have been performed on a representative slide of each lot of TMA product.

Tumor Microarray Map

Slide Label	TMA-HP-LI-063			triplicate								
	A	B	C	D	E	F	G	H	I	J	K	L
1	LI0348-P6 Liver cancer	LI0612-P12 Liver cancer	LI0752-P9 Liver cancer	LI0801-P6 Liver cancer	LI1035-P8 Liver cancer	LI1037-P11 Liver cancer	LI1057-P5 Liver cancer	LI1058-P7 Liver cancer	LI1068-P10 Liver cancer	LI1069-P6 Liver cancer	LI1078-P8 Liver cancer	LI1081-P5 Liver cancer
2	LI1097-P7 Liver cancer	LI1099-P8 Liver cancer	LI1549-P7 Liver cancer	LI1646-P3 Liver cancer	LI2658-P4 Liver cancer	LI3243-P5 Liver cancer	LI6281-P6 Liver cancer	LI6293-P10 Liver cancer	LI6306-P9 Liver cancer	LI6601-P6 Liver cancer	LI6610-P10 Liver cancer	LI6612-P8 Liver cancer
3	LI6616-P9 Liver cancer	LI6659-P8 Liver cancer	LI6669-P12 Liver cancer	LI6671-P9 Liver cancer	LI6677-P10 Liver cancer	LI6683-P7 Liver cancer	LI6722-P8 Liver cancer	LI9318-P2 Liver cancer	LI0348-P6 Liver cancer	LI0612-P12 Liver cancer	LI0752-P9 Liver cancer	LI0801-P6 Liver cancer
4	LI1035-P8 Liver cancer	LI1037-P11 Liver cancer	LI1057-P5 Liver cancer	LI1058-P7 Liver cancer	LI1068-P10 Liver cancer	LI1069-P6 Liver cancer	LI1078-P8 Liver cancer	LI1081-P5 Liver cancer	LI1097-P7 Liver cancer	LI1099-P8 Liver cancer	LI1549-P7 Liver cancer	LI1646-P3 Liver cancer
5	LI2658-P4 Liver cancer	LI3243-P5 Liver cancer	LI6281-P6 Liver cancer	LI6293-P10 Liver cancer	LI6306-P9 Liver cancer	LI6601-P6 Liver cancer	LI6610-P10 Liver cancer	LI6612-P8 Liver cancer	LI6616-P9 Liver cancer	LI6659-P8 Liver cancer	LI6669-P12 Liver cancer	LI6671-P9 Liver cancer
6	LI6677-P10 Liver cancer	LI6683-P7 Liver cancer	LI6722-P8 Liver cancer	LI9318-P2 Liver cancer	LI0348-P6 Liver cancer	LI0612-P12 Liver cancer	LI0752-0P9 Liver cancer	LI0801-P6 Liver cancer	LI1035-P8 Liver cancer	LI1037-P11 Liver cancer	LI1057-P5 Liver cancer	LI1058-P7 Liver cancer
7	LI1068-P10 Liver cancer	LI1069-P6 Liver cancer	LI1078-P8 Liver cancer	LI1081-P5 Liver cancer	LI1097-P7 Liver cancer	LI1099-P8 Liver cancer	LI1549-P7 Liver cancer	LI1646-P3 Liver cancer	LI2658-P4 Liver cancer	LI3243-P5 Liver cancer	LI6281-P6 Liver cancer	LI6293-P10 Liver cancer
8	LI6306-P9 Liver cancer	LI6601-P6 Liver cancer	LI6610-P10 Liver cancer	LI6612-P8 Liver cancer	LI6616-P9 Liver cancer	LI6659-P8 Liver cancer	LI6669-P12 Liver cancer	LI6671-P9 Liver cancer	LI6677-P10 Liver cancer	LI6683-P7 Liver cancer	LI6722-P8 Liver cancer	LI9318-P2 Liver cancer
												Mark point

Map of liver cancer tumor microarray slide

Please follow standard IHC staining procedure.

Visit <https://www.crownbio.com/databases/hubase> to view detailed information for each model. Access to HuBase™ requires online registration via the link provided.