



TMA048

Lung Cancer Tumor Microarrays derived from PDX for H&E and IHC staining, formalin-fixed and embedded in paraffin.

Tumor Microarray Profile

Catalog #: TMA-HP-LU-048

Model Type: PDX

Host Species: Mouse

Models On TMA: 32

Cores/Model: 3

Cores On TMA: 96

Core Size: 1 mm diameter

Fixation Method: Formalin fixed, paraffin embedded

Applications: H&E and IHC staining

Storage Conditions: 4°C Stable for 12 months from date of purchase when stored as recommended

Slide Preparation: Slides are sealed with wax; heat slide at 60° for 1 hour before de-waxing

QC: Each TMA slide has >90% core occupancy. H&E and Ki67 IHC staining have been performed on a representative slide of each lot of TMA product.

Tumor Microarray Map

Slide Label		TMA-HP-LU-048	triplicate										
1	LU11584-P3 Lung cancer	LU11592-P5 Lung cancer	LU11609-P2 Lung cancer	LU11672-P4 Lung cancer	LU11692-P4 Lung cancer	LU11709-P5 Lung cancer	LU11736-P2 Lung cancer	LU11841-P3 Lung cancer	LU11870-P3 Lung cancer	LU11876-P6 Lung cancer	LU11894-P7 Lung cancer	LU11942-P4 Lung cancer	
2	LU11951-P5 Lung cancer	LU13944-P5 Lung cancer	LU13953-P9 Lung cancer	LU5146-P5 Lung cancer	LU5147-P4 Lung cancer	LU5148-P6 Lung cancer	LU5150-P4 Lung cancer	LU5155-P6 Lung cancer	LU5157-P9 Lung cancer	LU5160-P4 Lung cancer	LU5169-P2 Lung cancer	LU5174-P7 Lung cancer	
3	LU5194-P3 Lung cancer	LU5212-P2 Lung cancer	LU5229-P8 Lung cancer	LU5237-P6 Lung cancer	LU5245-P6 Lung cancer	LU5247-P5 Lung cancer	LU5355-P7 Lung cancer	LU5371-P4 Lung cancer	LU11584-P3 Lung cancer	LU11592-P5 Lung cancer	LU11609-P2 Lung cancer	LU11672-P4 Lung cancer	
4	LU11692-P4 Lung cancer	LU11709-P5 Lung cancer	LU11736-P2 Lung cancer	LU11841-P3 Lung cancer	LU11870-P3 Lung cancer	LU11876-P6 Lung cancer	LU11894-P7 Lung cancer	LU11942-P4 Lung cancer	LU11951-P5 Lung cancer	LU13944-P5 Lung cancer	LU13953-P9 Lung cancer	LU5146-P5 Lung cancer	
5	LU5147-P4 Lung cancer	LU5148-P6 Lung cancer	LU5150-P4 Lung cancer	LU5155-P6 Lung cancer	LU5157-P9 Lung cancer	LU5160-P4 Lung cancer	LU5169-P2 Lung cancer	LU5174-P7 Lung cancer	LU5194-P3 Lung cancer	LU5212-P2 Lung cancer	LU5229-P8 Lung cancer	LU5237-P6 Lung cancer	
6	LU5245-P6 Lung cancer	LU5247-P5 Lung cancer	LU5355-P7 Lung cancer	LU5371-P4 Lung cancer	LU11584-P3 Lung cancer	LU11592-P5 Lung cancer	LU11609-P2 Lung cancer	LU11672-P4 Lung cancer	LU11692-P4 Lung cancer	LU11709-P5 Lung cancer	LU11736-P2 Lung cancer	LU11841-P3 Lung cancer	
7	LU11870-P3 Lung cancer	LU11876-P6 Lung cancer	LU11894-P7 Lung cancer	LU11942-P4 Lung cancer	LU11951-P5 Lung cancer	LU13944-P5 Lung cancer	LU13953-P9 Lung cancer	LU5146-P5 Lung cancer	LU5147-P4 Lung cancer	LU5148-P6 Lung cancer	LU5150-P4 Lung cancer	LU5155-P6 Lung cancer	
8	LU5157-P9 Lung cancer	LU5160-P4 Lung cancer	LU5169-P2 Lung cancer	LU5174-P7 Lung cancer	LU5194-P3 Lung cancer	LU5212-P2 Lung cancer	LU5229-P8 Lung cancer	LU5237-P6 Lung cancer	LU5245-P6 Lung cancer	LU5247-P5 Lung cancer	LU5355-P7 Lung cancer	LU5371-P4 Lung cancer	
													Mark point

Map of lung cancer tumor microarray

Please follow standard IHC staining procedure.

Visit <https://www.crownbio.com/databases/hubase> to view detailed information for each model. Access to HuBase™ requires online registration via the link provided.