



TMA058

Lung Cancer Microarrays derived from PDX for H&E and IHC staining, formalin-fixed and embedded in paraffin.

Tumor Microarray Profile

Catalog #: TMA-HP-LU-058

Model Type: PDX

Host Species: Mouse

Models On TMA: 32

Cores/Model: 3

Cores On TMA: 96

Core Size: 1 mm diameter

Fixation Method: Formalin fixed, paraffin embedded

Applications: H&E and IHC staining

Storage Conditions: 4°C Stable for 12 months from date of purchase when stored as recommended

Slide Preparation: Slides are sealed with wax; heat slide at 60° for 1 hour before de-waxing

QC: Each TMA slide has >90% core occupancy. H&E and Ki67 IHC staining have been performed on a representative slide of each lot of TMA product.

Tumor Microarray Map

Slide Label		TMA-HP-LU-058		triplicate									
	A	B	C	D	E	F	G	H	I	J	K	L	
1	LU0350-P4 Lung cancer	LU0755-P7 Lung cancer	LU0884-P8 Lung cancer	LU11622-P4 Lung cancer	LU11624-P8 Lung cancer	LU11784-P6 Lung cancer	LU11855-P5 Lung cancer	LU1656-P4 Lung cancer	LU1776-P9 Lung cancer	LU1828-P5 Lung cancer	LU3133-P3 Lung cancer	LU5149-P11 Lung cancer	
2	LU5173-P11 Lung cancer	LU5251-P9 Lung cancer	LU5351-P8 Lung cancer	LU5408-P6 Lung cancer	LU6402-P26 Lung cancer	LU6411-P19 Lung cancer	LU6423-P11 Lung cancer	LU6431-P16 Lung cancer	LU6434-P9 Lung cancer	LU6436-P9 Lung cancer	LU6437-P11 Lung cancer	LU6439-P7 Lung cancer	
3	LU6471-P10 Lung cancer	LU6481-P14 Lung cancer	LU6484-P12 Lung cancer	LU6800-P9 Lung cancer	LU6913-P9 Lung cancer	LU6917-P10 Lung cancer	LU6920-P9 Lung cancer	LU9902-P3 Lung cancer	LU0350-P4 Lung cancer	LU0755-P7 Lung cancer	LU0884-P8 Lung cancer	LU11622-P4 Lung cancer	
4	LU11624-P8 Lung cancer	LU11784-P6 Lung cancer	LU11855-P5 Lung cancer	LU1656-P4 Lung cancer	LU1776-P9 Lung cancer	LU1828-P5 Lung cancer	LU3133-P3 Lung cancer	LU5149-P11 Lung cancer	LU5173-P11 Lung cancer	LU5251-P9 Lung cancer	LU5351-P8 Lung cancer	LU5408-P6 Lung cancer	
5	LU6402-P26 Lung cancer	LU6411-P19 Lung cancer	LU6423-P11 Lung cancer	LU6431-P16 Lung cancer	LU6434-P9 Lung cancer	LU6436-P9 Lung cancer	LU6437-P11 Lung cancer	LU6439-P7 Lung cancer	LU6471-P10 Lung cancer	LU6481-P14 Lung cancer	LU6484-P12 Lung cancer	LU6800-P9 Lung cancer	
6	LU6913-P9 Lung cancer	LU6917-P10 Lung cancer	LU6920-P9 Lung cancer	LU9902-P3 Lung cancer	LU0350-P4 Lung cancer	LU0755-P7 Lung cancer	LU0884-P8 Lung cancer	LU11622-P4 Lung cancer	LU11624-P8 Lung cancer	LU11784-P6 Lung cancer	LU11855-P5 Lung cancer	LU1656-P4 Lung cancer	
7	LU1776-P9 Lung cancer	LU1828-P5 Lung cancer	LU3133-P3 Lung cancer	LU5149-P11 Lung cancer	LU5173-P11 Lung cancer	LU5251-P9 Lung cancer	LU5351-P8 Lung cancer	LU5408-P6 Lung cancer	LU6402-P26 Lung cancer	LU6411-P19 Lung cancer	LU6423-P11 Lung cancer	LU6431-P16 Lung cancer	
8	LU6434-P9 Lung cancer	LU6436-P9 Lung cancer	LU6437-P11 Lung cancer	LU6439-P7 Lung cancer	LU6471-P10 Lung cancer	LU6481-P14 Lung cancer	LU6484-P12 Lung cancer	LU6800-P9 Lung cancer	LU6913-P9 Lung cancer	LU6917-P10 Lung cancer	LU6920-P9 Lung cancer	LU9902-P3 Lung cancer	
												Mark point	

Map of lung cancer tumor microarray

Please follow standard IHC staining procedure.

Visit <https://www.crownbio.com/databases> to view detailed information for each model. Access to the databases requires online registration via the link provided.