



TMA059

Lung Cancer Microarrays derived from PDX for H&E and IHC staining, formalin-fixed and embedded in paraffin.

Tumor Microarray Profile

Catalog #: TMA-HP-LU-059

Model Type: PDX

Host Species: Mouse

Models On TMA: 32

Cores/Model: 3

Cores On TMA: 96

Core Size: 1 mm diameter

Fixation Method: Formalin fixed, paraffin embedded

Applications: H&E and IHC staining

Storage Conditions: 4°C Stable for 12 months from date of purchase when stored as recommended

Slide Preparation: Slides are sealed with wax; heat slide at 60° for 1 hour before de-waxing

QC: Each TMA slide has >90% core occupancy. H&E and Ki67 IHC staining have been performed on a representative slide of each lot of TMA product.

Tumor Microarray Map

Slide Label		TMA-HP-LU-059		triplicate									
	A	B	C	D	E	F	G	H	I	J	K	L	
1	LU11561-P5 Lung cancer	LU11568-P8 Lung cancer	LU11606-P7 Lung cancer	LU11634-P7 Lung cancer	LU11644-P7 Lung cancer	LU11956-P5 Lung cancer	LU1235-P5 Lung cancer	LU1380-P10 Lung cancer	LU2512-P7 Lung cancer	LU2535-P9 Lung cancer	LU2553-P9 Lung cancer	LU3057-P10 Lung cancer	
2	LU5151-P11 Lung cancer	LU5159-P8 Lung cancer	LU5176-P8 Lung cancer	LU5178-P10 Lung cancer	LU5187-P7 Lung cancer	LU5221-P6 Lung cancer	LU5268-P6 Lung cancer	LU5350-P9 Lung cancer	LU5407-P6 Lung cancer	LU6412-P13 Lung cancer	LU6424-P15 Lung cancer	LU6425-P10 Lung cancer	
3	LU6442-P8 Lung cancer	LU6802-P9 Lung cancer	LU6912-P8 Lung cancer	LU6918-P7 Lung cancer	LU9350-P7 Lung cancer	LU9543-P5 Lung cancer	LU9575-P4 Lung cancer	LU9661-P3 Lung cancer	LU11561-P5 Lung cancer	LU11568-P8 Lung cancer	LU11606-P7 Lung cancer	LU11634-P7 Lung cancer	
4	LU11644-P7 Lung cancer	LU11956-P5 Lung cancer	LU1235-P5 Lung cancer	LU1380-P10 Lung cancer	LU2512-P7 Lung cancer	LU2535-P9 Lung cancer	LU2553-P9 Lung cancer	LU3057-P10 Lung cancer	LU5151-P11 Lung cancer	LU5159-P8 Lung cancer	LU5176-P8 Lung cancer	LU5178-P10 Lung cancer	
5	LU5187-P7 Lung cancer	LU5221-P6 Lung cancer	LU5268-P6 Lung cancer	LU5350-P9 Lung cancer	LU5407-P6 Lung cancer	LU6412-P13 Lung cancer	LU6424-P15 Lung cancer	LU6425-P10 Lung cancer	LU6442-P8 Lung cancer	LU6802-P9 Lung cancer	LU6912-P8 Lung cancer	LU6918-P7 Lung cancer	
6	LU9350-P7 Lung cancer	LU9543-P5 Lung cancer	LU9575-P4 Lung cancer	LU9661-P3 Lung cancer	LU11561-P5 Lung cancer	LU11568-P8 Lung cancer	LU11606-P7 Lung cancer	LU11634-P7 Lung cancer	LU11644-P7 Lung cancer	LU11956-P5 Lung cancer	LU1235-P5 Lung cancer	LU1380-P10 Lung cancer	
7	LU2512-P7 Lung cancer	LU2535-P9 Lung cancer	LU2553-P9 Lung cancer	LU3057-P10 Lung cancer	LU5151-P11 Lung cancer	LU5159-P8 Lung cancer	LU5176-P8 Lung cancer	LU5178-P10 Lung cancer	LU5187-P7 Lung cancer	LU5221-P6 Lung cancer	LU5268-P6 Lung cancer	LU5350-P9 Lung cancer	
8	LU5407-P6 Lung cancer	LU6412-P13 Lung cancer	LU6424-P15 Lung cancer	LU6425-P10 Lung cancer	LU6442-P8 Lung cancer	LU6802-P9 Lung cancer	LU6912-P8 Lung cancer	LU6918-P7 Lung cancer	LU9350-P7 Lung cancer	LU9543-P5 Lung cancer	LU9575-P4 Lung cancer	LU9661-P3 Lung cancer	
												Mark point	

Map of lung cancer tumor microarray

Please follow standard IHC staining procedure.

Visit <https://www.crownbio.com/databases> to view detailed information for each model. Access to the databases requires online registration via the link provided.