



TMA067

Lung Cancer Microarrays derived from PDX for H&E and IHC staining, formalin-fixed and embedded in paraffin.

Tumor Microarray Profile

Catalog #: TMA-HP-LU-067

Model Type: PDX

Host Species: Mouse

Models On TMA: 32

Cores/Model: 3

Cores On TMA: 96

Core Size: 1 mm diameter

Fixation Method: Formalin fixed, paraffin embedded

Applications: H&E and IHC staining

Storage Conditions: 4°C Stable for 12 months from date of purchase when stored as recommended

Slide Preparation: Slides are sealed with wax; heat slide at 60° for 1 hour before de-waxing

QC: Each TMA slide has >90% core occupancy. H&E and Ki67 IHC staining have been performed on a representative slide of each lot of TMA product.

Tumor Microarray Map

Slide Label		TMA-HP-LU-067	triplicate										
	A	B	C	D	E	F	G	H	I	J	K	L	
1	LU0748-P4 Lung cancer	LU0861-P10 Lung cancer	LU0876-P11 Lung cancer	LU0894-P6 Lung cancer	LU1116-P8 Lung cancer	LU1144-P7 Lung cancer	LU11693-P7 Lung cancer	LU11856-P8 Lung cancer	LU11885-P7 Lung cancer	LU1429-P9 Lung cancer	LU1452-P9 Lung cancer	LU1690-P7 Lung cancer	
2	LU20004-P2 Lung cancer	LU20123-P5 Lung cancer	LU2086-P10 Lung cancer	LU2555-P8 Lung cancer	LU2569-P8 Lung cancer	LU3073-P7 Lung cancer	LU5139-P10 Lung cancer	LU5185-P12 Lung cancer	LU5209-P8 Lung cancer	LU6266-P7 Lung cancer	LU6307-P10 Lung cancer	LU6407-P18 Lung cancer	
3	LU6414-P10 Lung cancer	LU6419-P12 Lung cancer	LU6907-P9 Lung cancer	LU6914-P11 Lung cancer	LU6919-P9 Lung cancer	LU6922-P11 Lung cancer	LU9536-P3 Lung cancer	LU9906-P5 Lung cancer	LU0748-P4 Lung cancer	LU0861-P10 Lung cancer	LU0876-P11 Lung cancer	LU0894-P6 Lung cancer	
4	LU1116-P8 Lung cancer	LU1144-P7 Lung cancer	LU11693-P7 Lung cancer	LU11856-P8 Lung cancer	LU11885-P7 Lung cancer	LU1429-P9 Lung cancer	LU1452-P9 Lung cancer	LU1690-P7 Lung cancer	LU20004-P2 Lung cancer	LU20123-P5 Lung cancer	LU2086-P10 Lung cancer	LU2555-P8 Lung cancer	
5	LU2569-P8 Lung cancer	LU3073-P7 Lung cancer	LU5139-P10 Lung cancer	LU5185-P12 Lung cancer	LU5209-P8 Lung cancer	LU6266-P7 Lung cancer	LU6307-P10 Lung cancer	LU6407-P18 Lung cancer	LU6414-P10 Lung cancer	LU6419-P12 Lung cancer	LU6907-P9 Lung cancer	LU6914-P11 Lung cancer	
6	LU6919-P9 Lung cancer	LU6922-P11 Lung cancer	LU9536-P3 Lung cancer	LU9906-P5 Lung cancer	LU0748-P4 Lung cancer	LU0861-P10 Lung cancer	LU0876-P11 Lung cancer	LU0894-P6 Lung cancer	LU1116-P8 Lung cancer	LU1144-P7 Lung cancer	LU11693-P7 Lung cancer	LU11856-P8 Lung cancer	
7	LU11885-P7 Lung cancer	LU1429-P9 Lung cancer	LU1452-P9 Lung cancer	LU1690-P7 Lung cancer	LU20004-P2 Lung cancer	LU20123-P5 Lung cancer	LU2086-P10 Lung cancer	LU2555-P8 Lung cancer	LU2569-P8 Lung cancer	LU3073-P7 Lung cancer	LU5139-P10 Lung cancer	LU5185-P12 Lung cancer	
8	LU5209-P8 Lung cancer	LU6266-P7 Lung cancer	LU6307-P10 Lung cancer	LU6407-P18 Lung cancer	LU6414-P10 Lung cancer	LU6419-P12 Lung cancer	LU6907-P9 Lung cancer	LU6914-P11 Lung cancer	LU6919-P9 Lung cancer	LU6922-P11 Lung cancer	LU9536-P3 Lung cancer	LU9906-P5 Lung cancer	
												Mark point	

Map of lung cancer tumor microarray

Please follow standard IHC staining procedure.

Visit <https://www.crownbio.com/databases> to view detailed information for each model. Access to the databases requires online registration via the link provided.