



TMA029

This ovarian cancer tumor microarray is derived from patient derived xenograft (PDX) tumor samples. Each tumor sample is formalin fixed, paraffin embedded, and spotted in triplicate on the TMA.

Tumor Microarray Profile

Catalog #: TMA-HP-OV-029

Model Type: PDX

Host Species: Mouse

Models On TMA: 28

Cores/Model: 3

Cores On TMA: 84

Core Size: 1 mm diameter

Fixation Method: Formalin fixed, paraffin embedded

Applications: H&E and IHC staining

Storage Conditions: 4°C Stable for 12 months from date of purchase when stored as recommended

Slide Preparation: Slides are sealed with wax; heat slide at 60° for 1 hour before de-waxing

QC: Each TMA slide has >90% core occupancy. H&E and Ki67 IHC staining have been performed on a representative slide of each lot of TMA product.

Tumor Microarray Map

Slide Label	TMA-HP-OV-029												
	triplicate												
	A	B	C	D	E	F	G	H	I	J	K	L	
	1	OV0243-P5 Ovarian cancer	OV0250-P11 Ovarian cancer	OV0273-P9 Ovarian cancer	OV0276-P5 Ovarian cancer	OV0589-P2 Ovarian cancer	OV1200-P4 Ovarian cancer	OV1272-P13 Ovarian cancer	OV1286-P7 Ovarian cancer	OV1327-P7 Ovarian cancer	OV1385-P4 Ovarian cancer	OV1579-P6 Ovarian cancer	OV1658-P9 Ovarian cancer
	2	OV2015-P4 Ovarian cancer	OV2017-P10 Ovarian cancer	OV2018-P5 Ovarian cancer	OV2026-P5 Ovarian cancer	OV2027-P7 Ovarian cancer	OV2029-P12 Ovarian cancer	OV2034-P11 Ovarian cancer	OV2423-P10 Ovarian cancer	OV2546-P7 Ovarian cancer	OV3355-P3 Ovarian cancer	OV3395-P5 Ovarian cancer	OV3414-P3 Ovarian cancer
	3	OV3425-P2 Ovarian cancer	OV3596-P6 Ovarian cancer	OV3756-P3 Ovarian cancer	unknown	OV6308-P4 Ovarian cancer	OV0243-P5 Ovarian cancer	OV0250-P11 Ovarian cancer	OV0273-P9 Ovarian cancer	OV0276-P5 Ovarian cancer	OV0589-P2 Ovarian cancer	OV1200-P4 Ovarian cancer	OV1272-P13 Ovarian cancer
	4	OV1286-P7 Ovarian cancer	OV1327-P7 Ovarian cancer	OV1385-P4 Ovarian cancer	OV1579-P6 Ovarian cancer	OV1658-P9 Ovarian cancer	OV2015-P4 Ovarian cancer	OV2017-P10 Ovarian cancer	OV2018-P5 Ovarian cancer	OV2026-P5 Ovarian cancer	OV2027-P7 Ovarian cancer	OV2029-P12 Ovarian cancer	OV2034-P11 Ovarian cancer
	5	OV2423-P10 Ovarian cancer	OV2546-P7 Ovarian cancer	OV3355-P3 Ovarian cancer	OV3395-P5 Ovarian cancer	OV3414-P3 Ovarian cancer	OV3425-P2 Ovarian cancer	OV3596-P6 Ovarian cancer	OV3756-P3 Ovarian cancer	unknown	OV6308-P4 Ovarian cancer	OV0243-P5 Ovarian cancer	OV0250-P11 Ovarian cancer
	6	OV0273-P9 Ovarian cancer	OV0276-P5 Ovarian cancer	OV0589-P2 Ovarian cancer	OV1200-P4 Ovarian cancer	OV1272-P13 Ovarian cancer	OV1286-P7 Ovarian cancer	OV1327-P7 Ovarian cancer	OV1385-P4 Ovarian cancer	OV1579-P6 Ovarian cancer	OV1658-P9 Ovarian cancer	OV2015-P4 Ovarian cancer	OV2017-P10 Ovarian cancer
	7	OV2018-P5 Ovarian cancer	OV2026-P5 Ovarian cancer	OV2027-P7 Ovarian cancer	OV2029-P12 Ovarian cancer	OV2034-P11 Ovarian cancer	OV2423-P10 Ovarian cancer	OV2546-P7 Ovarian cancer	OV3355-P3 Ovarian cancer	OV3395-P5 Ovarian cancer	OV3414-P3 Ovarian cancer	OV3425-P2 Ovarian cancer	OV3596-P6 Ovarian cancer
8	OV3756-P3 Ovarian cancer	unknown	OV6308-P4 Ovarian cancer										
												Mark point	

Map of breast cancer tumor microarray

Please follow standard IHC staining procedure.

Visit <https://www.crownbio.com/databases/hubase> to view detailed information for each model. Access to HuBase™ requires online registration via the link provided.