

TMA060

Ovarian Cancer Microarrays derived from PDX for H&E and IHC staining, formalin-fixed and embedded in paraffin.

Tumor Microarray Profile

Catalog #: TMA-HP-OV-060

Model Type: PDX

Host Species: Mouse

Models On TMA: 32

Cores/Model: 3

Cores On TMA:96

Core Size: 1 mm diameter

Fixation Method: Formalin fixed, paraffin embedded

Applications: H&E and IHC staining

Storage Conditions: 4°C Stable for 12 months from date of purchase when stored as recommended

Slide Preparation: Slides are sealed with wax; heat slide at 60° for 1 hour before de-waxing

QC: Each TMA slide has >90% core occupancy. H&E and Ki67 IHC staining have been performed on a representative slide of each lot of TMA product.

Tumor Microarray Map

Slide		TMA-HP-OV-0	60	triplicate									
Label		Α	В	С	D	E	F	G	Н	- 1	J	K	L
	1	OV13951-P8 Ovarian cancer	OV14108-P5 Ovarian cancer	OV14930-P5 Ovarian cancer	OV14948-P6 Ovarian cancer	OV15209-P11 Ovarian cancer	OV15387-P6 Ovarian cancer	OV15631-P10 Ovarian cancer	OV15727-P4 Ovarian cancer	OV3708-P7 Ovarian cancer	OV5287-P6 Ovarian cancer	OV5294-P8 Ovarian cancer	OV5296-P7 Ovarian cancer
	2	OV5297-P2 Ovarian cancer	OV5304-P6 Ovarian cancer	OV5390-P12 Ovarian cancer	OV9401-P7 Ovarian cancer	OV9402-P6 Ovarian cancer	OV9403-P3 Ovarian cancer	OV9404-P5 Ovarian cancer	OV9412-P10 Ovarian cancer	OV9419-P10 Ovarian cancer	OV9420-P7 Ovarian cancer	OV9421-P10 Ovarian cancer	OV9424-P7 Ovarian cancer
	3	OV9427-P7 Ovarian cancer	OV9428-P12 Ovarian cancer	OV9497-P8 Ovarian cancer	OV9502-P12 Ovarian cancer	OV9522-P10 Ovarian cancer	OV9540-P7 Ovarian cancer	OV9579-P4 Ovarian cancer	OV9580-P3 Ovarian cancer	OV13951-P8 Ovarian cancer	OV14108-P5 Ovarian cancer	OV14930-P5 Ovarian cancer	OV14948-P6 Ovarian cancer
	4	OV15209-P11 Ovarian cancer	OV15387-P6 Ovarian cancer	OV15631-P10 Ovarian cancer	OV15727-P4 Ovarian cancer	OV3708-P7 Ovarian cancer	OV5287-P6 Ovarian cancer	OV5294-P8 Ovarian cancer	OV5296-P7 Ovarian cancer	OV5297-P2 Ovarian cancer	OV5304-P6 Ovarian cancer	OV5390-P12 Ovarian cancer	OV9401-P7 Ovarian cancer
	5	OV9402-P6 Ovarian cancer	OV9403-P3 Ovarian cancer	OV9404-P5 Ovarian cancer	OV9412-P10 Ovarian cancer	OV9419-P10 Ovarian cancer	OV9420-P7 Ovarian cancer	OV9421-P10 Ovarian cancer	OV9424-P7 Ovarian cancer	OV9427-P7 Ovarian cancer	OV9428-P12 Ovarian cancer	OV9497-P8 Ovarian cancer	OV9502-P12 Ovarian cancer
	6	OV9522-P10 Ovarian cancer	OV9540-P7 Ovarian cancer	OV9579-P4 Ovarian cancer	OV9580-P3 Ovarian cancer	OV13951-P8 Ovarian cancer	OV14108-P5 Ovarian cancer	OV14930-P5 Ovarian cancer	OV14948-P6 Ovarian cancer	OV15209-P11 Ovarian cancer	OV15387-P6 Ovarian cancer	OV15631-P10 Ovarian cancer	OV15727-P4 Ovarian cancer
	7	OV3708-P7 Ovarian cancer	OV5287-P6 Ovarian cancer	OV5294-P8 Ovarian cancer	OV5296-P7 Ovarian cancer	OV5297-P2 Ovarian cancer	OV5304-P6 Ovarian cancer	OV5390-P12 Ovarian cancer	OV9401-P7 Ovarian cancer	OV9402-P6 Ovarian cancer	OV9403-P3 Ovarian cancer	OV9404-P5 Ovarian cancer	OV9412-P10 Ovarian cancer
	8	OV9419-P10 Ovarian cancer	OV9420-P7 Ovarian cancer	OV9421-P10 Ovarian cancer	OV9424-P7 Ovarian cancer	OV9427-P7 Ovarian cancer	OV9428-P12 Ovarian cancer	OV9497-P8 Ovarian cancer	OV9502-P12 Ovarian cancer	OV9522-P10 Ovarian cancer	OV9540-P7 Ovarian cancer	OV9579-P4 Ovarian cancer	OV9580-P3 Ovarian cancer
													Mark point

Map of ovarian cancer tumor microarray

Please follow standard IHC staining procedure.

Visit https://www.crownbio.com/databases to view detailed information for each model. Access to the databases requires online registration via the link provided.