



TMA018

This pancreatic cancer tumor microarray is derived from patient derived xenograft (PDX) tumor samples. Each tumor sample is formalin fixed, paraffin embedded, and spotted in triplicate on the TMA.

Tumor Microarray Profile

Catalog #: TMA-HP-PA-018

Model Type: PDX

Host Species: Mouse

Models On TMA: 31

Cores/Model: 3

Cores On TMA: 93

Core Size: 1 mm diameter

Fixation Method: Formalin fixed, paraffin embedded

Applications: H&E and IHC staining

Storage Conditions: 4°C Stable for 12 months from date of purchase when stored as recommended

Slide Preparation: Slides are sealed with wax; heat slide at 60° for 1 hour before de-waxing

QC: Each TMA slide has >90% core occupancy. H&E and Ki67 IHC staining have been performed on a representative slide of each lot of TMA product.

Tumor Microarray Map

Slide Label	TMA-HP-PA-018											
	A	B	C	D	E	F	G	H	I	J	K	L
1	PA6285-P13 Pancreatic cancer	PA6259-P4 Pancreatic cancer	PA6236-P3 Pancreatic cancer	PA3065-P5 Pancreatic cancer	PA1443-P7 Pancreatic cancer	unknown-NA unknown	PA1383-P8 Pancreatic cancer	PA1367-P3 Pancreatic cancer	PA1360-P6 Pancreatic cancer	PA1357-P4 Pancreatic cancer	PA1338-P2 Pancreatic cancer	PA1337-P6 Pancreatic cancer
2	PA1332-P9 Pancreatic cancer	PA1316-P4 Pancreatic cancer	PA1306-P5 Pancreatic cancer	PA1301-P4 Pancreatic cancer	PA1266-P6 Pancreatic cancer	PA1252-P3 Pancreatic cancer	PA1233-P5 Pancreatic cancer	PA1280-P3 Pancreatic cancer	PA1229-P7 Pancreatic cancer	PA1222-P4 Pancreatic cancer	PA1221-P6 Pancreatic cancer	PA1214-P4 Pancreatic cancer
3	PA1198-P7 Pancreatic cancer	PA1194-P10 Pancreatic cancer	PA1189-P5 Pancreatic cancer	PA1178-P10 Pancreatic cancer	PA1170-P6 Pancreatic cancer	PA1168-P7 Pancreatic cancer	PA0787-P4 Pancreatic cancer	PA0372-P5 Pancreatic cancer	PA6285-P13 Pancreatic cancer	PA6259-P4 Pancreatic cancer	PA6236-P3 Pancreatic cancer	PA3065-P5 Pancreatic cancer
4	PA1443-P7 Pancreatic cancer	unknown-NA unknown	PA1383-P8 Pancreatic cancer	PA1367-P3 Pancreatic cancer	PA1360-P6 Pancreatic cancer	PA1357-P4 Pancreatic cancer	PA1338-P2 Pancreatic cancer	PA1337-P6 Pancreatic cancer	PA1332-P9 Pancreatic cancer	PA1316-P4 Pancreatic cancer	PA1306-P5 Pancreatic cancer	PA1301-P4 Pancreatic cancer
5	PA1266-P6 Pancreatic cancer	PA1252-P3 Pancreatic cancer	PA1233-P5 Pancreatic cancer	PA1280-P3 Pancreatic cancer	PA1229-P7 Pancreatic cancer	PA1222-P4 Pancreatic cancer	PA1221-P6 Pancreatic cancer	PA1214-P4 Pancreatic cancer	PA1198-P7 Pancreatic cancer	PA1194-P10 Pancreatic cancer	PA1189-P5 Pancreatic cancer	PA1178-P10 Pancreatic cancer
6	PA1170-P6 Pancreatic cancer	PA1168-P7 Pancreatic cancer	PA0787-P4 Pancreatic cancer	PA0372-P5 Pancreatic cancer	PA6285-P13 Pancreatic cancer	PA6259-P4 Pancreatic cancer	PA6236-P3 Pancreatic cancer	PA3065-P5 Pancreatic cancer	PA1443-P7 Pancreatic cancer	unknown-NA unknown	PA1383-P8 Pancreatic cancer	PA1367-P3 Pancreatic cancer
7	PA1360-P6 Pancreatic cancer	PA1357-P4 Pancreatic cancer	PA1338-P2 Pancreatic cancer	PA1337-P6 Pancreatic cancer	PA1332-P9 Pancreatic cancer	PA1316-P4 Pancreatic cancer	PA1306-P5 Pancreatic cancer	PA1301-P4 Pancreatic cancer	PA1266-P6 Pancreatic cancer	PA1252-P3 Pancreatic cancer	PA1233-P5 Pancreatic cancer	PA1280-P3 Pancreatic cancer
8	PA1229-P7 Pancreatic cancer	PA1222-P4 Pancreatic cancer	PA1221-P6 Pancreatic cancer	PA1214-P4 Pancreatic cancer	PA1198-P7 Pancreatic cancer	PA1194-P10 Pancreatic cancer	PA1189-P5 Pancreatic cancer	PA1178-P10 Pancreatic cancer	PA1170-P6 Pancreatic cancer	PA1168-P7 Pancreatic cancer	PA0787-P4 Pancreatic cancer	PA0372-P5 Pancreatic cancer
												Mark point

Map of pancreatic cancer tumor microarray slide

Please follow standard IHC staining procedure.

Visit <https://www.crownbio.com/databases/hubase> to view detailed information for each model. Access to HuBase™ requires online registration via the link provided.