



TMA039

This pancreatic cancer tumor microarray is derived from patient derived xenograft (PDX) tumor samples. Each tumor sample is formalin fixed, paraffin embedded, and spotted in triplicate on the TMA.

Tumor Microarray Profile

Catalog #: TMA-HP-PA-039

Model Type: PDX

Host Species: Mouse

Models On TMA: 32

Cores/Model: 3

Cores On TMA: 96

Core Size: 1 mm diameter

Fixation Method: Formalin fixed, paraffin embedded

Applications: H&E and IHC staining

Storage Conditions: 4°C Stable for 12 months from date of purchase when stored as recommended

Slide Preparation: Slides are sealed with wax; heat slide at 60° for 1 hour before de-waxing

QC: Each TMA slide has >90% core occupancy. H&E and Ki67 IHC staining have been performed on a representative slide of each lot of TMA product.

Tumor Microarray Map

Slide Label	TMA-HP-PA-039													
	A	B	triplicate			C	D	E	F	G	H	I	J	K
1	PA9452-P2 Pancreatic cancer	PA9450-P2 Pancreatic cancer	PA6265-P10 Pancreatic cancer	PA6263-P3 Pancreatic cancer	PA6262-P8 Pancreatic cancer	PA6260-P4 Pancreatic cancer	PA6258-P7 Pancreatic cancer	PA6247-P5 Pancreatic cancer	PA6245-P5 Pancreatic cancer	PA6238-P12 Pancreatic cancer	PA6237-P5 Pancreatic cancer	PA6230-P10 Pancreatic cancer		
2	PA3142-P6 Pancreatic cancer	PA3139-P7 Pancreatic cancer	PA3137-P4 Pancreatic cancer	PA3126-P4 Pancreatic cancer	PA3060-P10 Pancreatic cancer	PA3029-P11 Pancreatic cancer	PA3013-P3 Pancreatic cancer	PA3008-P11 Pancreatic cancer	PA3004-P8 Pancreatic cancer	PA3003-P7 Pancreatic cancer	PA2637-P8 Pancreatic cancer	PA2566-P8 Pancreatic cancer		
3	PA1459-P9 Pancreatic cancer	PA1457-P6 Pancreatic cancer	PA1455-P9 Pancreatic cancer	PA1413-P10 Pancreatic cancer	PA1405-P9 Pancreatic cancer	PA1265-P5 Pancreatic cancer	PA0968-P3 Pancreatic cancer	PA0692-P9 Pancreatic cancer	PA9452-P2 Pancreatic cancer	PA9450-P2 Pancreatic cancer	PA6265-P10 Pancreatic cancer	PA6263-P3 Pancreatic cancer		
4	PA6262-P8 Pancreatic cancer	PA6260-P4 Pancreatic cancer	PA6258-P7 Pancreatic cancer	PA6247-P5 Pancreatic cancer	PA6245-P5 Pancreatic cancer	PA6238-P12 Pancreatic cancer	PA6237-P5 Pancreatic cancer	PA6230-P10 Pancreatic cancer	PA3142-P6 Pancreatic cancer	PA3139-P7 Pancreatic cancer	PA3137-P4 Pancreatic cancer	PA3126-P4 Pancreatic cancer		
5	PA3060-P10 Pancreatic cancer	PA3029-P11 Pancreatic cancer	PA3013-P3 Pancreatic cancer	PA3008-P11 Pancreatic cancer	PA3004-P8 Pancreatic cancer	PA3003-P7 Pancreatic cancer	PA2637-P8 Pancreatic cancer	PA2566-P8 Pancreatic cancer	PA1459-P9 Pancreatic cancer	PA1457-P6 Pancreatic cancer	PA1455-P9 Pancreatic cancer	PA1413-P10 Pancreatic cancer		
6	PA1405-P9 Pancreatic cancer	PA1265-P5 Pancreatic cancer	PA0968-P3 Pancreatic cancer	PA0692-P9 Pancreatic cancer	PA9452-P2 Pancreatic cancer	PA9450-P2 Pancreatic cancer	PA6265-P10 Pancreatic cancer	PA6263-P3 Pancreatic cancer	PA6262-P8 Pancreatic cancer	PA6260-P4 Pancreatic cancer	PA6258-P7 Pancreatic cancer	PA6247-P5 Pancreatic cancer		
7	PA6245-P5 Pancreatic cancer	PA6238-P12 Pancreatic cancer	PA6237-P5 Pancreatic cancer	PA6230-P10 Pancreatic cancer	PA3142-P6 Pancreatic cancer	PA3139-P7 Pancreatic cancer	PA3137-P4 Pancreatic cancer	PA3126-P4 Pancreatic cancer	PA3060-P10 Pancreatic cancer	PA3029-P11 Pancreatic cancer	PA3013-P3 Pancreatic cancer	PA3008-P11 Pancreatic cancer		
8	PA3004-P8 Pancreatic cancer	PA3003-P7 Pancreatic cancer	PA2637-P8 Pancreatic cancer	PA2566-P8 Pancreatic cancer	PA1459-P9 Pancreatic cancer	PA1457-P6 Pancreatic cancer	PA1455-P9 Pancreatic cancer	PA1413-P10 Pancreatic cancer	PA1405-P9 Pancreatic cancer	PA1265-P5 Pancreatic cancer	PA0968-P3 Pancreatic cancer	PA0692-P9 Pancreatic cancer		
												Mark point		

Map of pancreatic cancer tumor microarray slide

Please follow standard IHC staining procedure.

Visit <https://www.crownbio.com/databases/hubase> to view detailed information for each model. Access to HuBase™ requires online registration via the link provided.