



TMA06 1

Pancreatic Cancer Microarrays derived from PDX for H&E and IHC staining, formalin-fixed and embedded in paraffin.

Tumor Microarray Profile

Catalog #:TMA-HP-PA-061

Model Type:PDX

Host Species:Mouse

Models On TMA:32

Cores/Model:3

Cores On TMA:96

Core Size:1 mm diameter

Fixation Method:Formalin fixed, paraffin embedded

Applications:H&E and IHC staining

Storage Conditions:4°C Stable for 12 months from date of purchase when stored as recommended

Slide Preparation:Slides are sealed with wax; heat slide at 60° for 1 hour before de-waxing

QC:Each TMA slide has >90% core occupancy. H&E and Ki67 IHC staining have been performed on a representative slide of each lot of TMA product.

Tumor Microarray Map

Slide Label		TMA-HP-PA-061		triplicate									
		A	B	C	D	E	F	G	H	I	J	K	L
1	1	PA13026-P5 Pancreatic cancer	PA13033-P5 Pancreatic cancer	PA13034-P4 Pancreatic cancer	PA20006-P2 Pancreatic cancer	PA20066-P2 Pancreatic cancer	PA20067-P2 Pancreatic cancer	PA20068-P3 Pancreatic cancer	PA20069-P3 Pancreatic cancer	PA20070-P2 Pancreatic cancer	PA20072-P2 Pancreatic cancer	PA20073-P3 Pancreatic cancer	PA20075-P4 Pancreatic cancer
	2	PA20076-P2 Pancreatic cancer	PA20082-P3 Pancreatic cancer	PA20084-P4 Pancreatic cancer	PA20139-P6 Pancreatic cancer	PA20166-P7 Pancreatic cancer	PA3127-P7 Pancreatic cancer	PA3136-P11 Pancreatic cancer	PA3258-P3 Pancreatic cancer	PA5363-P2 Pancreatic cancer	PA5389-P4 Pancreatic cancer	PA5404-P3 Pancreatic cancer	PA5410-P2 Pancreatic cancer
	3	PA5411-P5 Pancreatic cancer	PA5412-P8 Pancreatic cancer	PA5415-P5 Pancreatic cancer	PA9435-P3 Pancreatic cancer	PA9449-P4 Pancreatic cancer	PA9451-P12 Pancreatic cancer	PA9453-P4 Pancreatic cancer	PA9573-P6 Pancreatic cancer	PA13026-P5 Pancreatic cancer	PA13033-P5 Pancreatic cancer	PA13034-P4 Pancreatic cancer	PA20006-P2 Pancreatic cancer
	4	PA20066-P2 Pancreatic cancer	PA20067-P2 Pancreatic cancer	PA20068-P3 Pancreatic cancer	PA20069-P3 Pancreatic cancer	PA20070-P2 Pancreatic cancer	PA20072-P2 Pancreatic cancer	PA20073-P3 Pancreatic cancer	PA20075-P4 Pancreatic cancer	PA20076-P2 Pancreatic cancer	PA20082-P3 Pancreatic cancer	PA20084-P4 Pancreatic cancer	PA20139-P6 Pancreatic cancer
	5	PA20166-P7 Pancreatic cancer	PA3127-P7 Pancreatic cancer	PA3136-P11 Pancreatic cancer	PA3258-P3 Pancreatic cancer	PA5363-P2 Pancreatic cancer	PA5389-P4 Pancreatic cancer	PA5404-P3 Pancreatic cancer	PA5410-P2 Pancreatic cancer	PA5411-P5 Pancreatic cancer	PA5412-P8 Pancreatic cancer	PA5415-P5 Pancreatic cancer	PA9435-P3 Pancreatic cancer
	6	PA9449-P4 Pancreatic cancer	PA9451-P12 Pancreatic cancer	PA9453-P4 Pancreatic cancer	PA9573-P6 Pancreatic cancer	PA13026-P5 Pancreatic cancer	PA13033-P5 Pancreatic cancer	PA13034-P4 Pancreatic cancer	PA20006-P2 Pancreatic cancer	PA20066-P2 Pancreatic cancer	PA20067-P2 Pancreatic cancer	PA20068-P3 Pancreatic cancer	PA20069-P3 Pancreatic cancer
	7	PA20070-P2 Pancreatic cancer	PA20072-P2 Pancreatic cancer	PA20073-P3 Pancreatic cancer	PA20075-P4 Pancreatic cancer	PA20076-P2 Pancreatic cancer	PA20082-P3 Pancreatic cancer	PA20084-P4 Pancreatic cancer	PA20139-P6 Pancreatic cancer	PA20166-P7 Pancreatic cancer	PA3127-P7 Pancreatic cancer	PA3136-P11 Pancreatic cancer	PA3258-P3 Pancreatic cancer
	8	PA5363-P2 Pancreatic cancer	PA5389-P4 Pancreatic cancer	PA5404-P3 Pancreatic cancer	PA5410-P2 Pancreatic cancer	PA5411-P5 Pancreatic cancer	PA5412-P8 Pancreatic cancer	PA5415-P5 Pancreatic cancer	PA9435-P3 Pancreatic cancer	PA9449-P4 Pancreatic cancer	PA9451-P12 Pancreatic cancer	PA9453-P4 Pancreatic cancer	PA9573-P6 Pancreatic cancer
													Mark point

Map of pancreatic cancer tumor microarray

Please follow standard IHC staining procedure.

Visit <https://www.crownbio.com/databases> to view detailed information for each model. Access to the databases requires online registration via the link provided.