



TMA062

Pancreatic Cancer Microarrays derived from PDX for H&E and IHC staining, formalin-fixed and embedded in paraffin.

Tumor Microarray Profile

Catalog #:TMA-HP-PA-062

Model Type:PDX

Host Species:Mouse

Models On TMA:32

Cores/Model:3

Cores On TMA:96

Core Size:1 mm diameter

Fixation Method:Formalin fixed, paraffin embedded

Applications:H&E and IHC staining

Storage Conditions:4°C Stable for 12 months from date of purchase when stored as recommended

Slide Preparation:Slides are sealed with wax; heat slide at 60° for 1 hour before de-waxing

QC:Each TMA slide has >90% core occupancy. H&E and Ki67 IHC staining have been performed on a representative slide of each lot of TMA product.

Tumor Microarray Map

Slide Label		TMA-HP-PA-062	triplicate										
	A	B	C	D	E	F	G	H	I	J	K	L	
1	PA0372-P5 Pancreatic cancer	PA0787-P5 Pancreatic cancer	PA1168-P4 Pancreatic cancer	PA1170-P2 Pancreatic cancer	PA1194-P2 Pancreatic cancer	PA1198-P2 Pancreatic cancer	PA1221-P3 Pancreatic cancer	PA1222-P3 Pancreatic cancer	PA1229-P2 Pancreatic cancer	PA1252-P2 Pancreatic cancer	PA1280-P3 Pancreatic cancer	PA1301-P4 Pancreatic cancer	
2	PA1332-P2 Pancreatic cancer	PA1337-P3 Pancreatic cancer	PA1338-P4 Pancreatic cancer	PA1360-P6 Pancreatic cancer	PA1383-P7 Pancreatic cancer	PA1443-P7 Pancreatic cancer	PA6236-P11 Pancreatic cancer	PA6285-P3 Pancreatic cancer	PA1290-P2 Pancreatic cancer	PA1322-P4 Pancreatic cancer	PA1390-P3 Pancreatic cancer	PA2021-P2 Pancreatic cancer	
3	PA2022-P5 Pancreatic cancer	PA2031-P8 Pancreatic cancer	PA2425-P5 Pancreatic cancer	PA3546-P3 Pancreatic cancer	PA6284-P4 Pancreatic cancer	PA2847-P12 Pancreatic cancer	PA2924-P4 Pancreatic cancer	PA3023-P6 Pancreatic cancer	PA0372-P5 Pancreatic cancer	PA0787-P5 Pancreatic cancer	PA1168-P4 Pancreatic cancer	PA1170-P2 Pancreatic cancer	
4	PA1194-P2 Pancreatic cancer	PA1198-P2 Pancreatic cancer	PA1221-P3 Pancreatic cancer	PA1222-P3 Pancreatic cancer	PA1229-P2 Pancreatic cancer	PA1252-P2 Pancreatic cancer	PA1280-P3 Pancreatic cancer	PA1301-P4 Pancreatic cancer	PA1332-P2 Pancreatic cancer	PA1337-P3 Pancreatic cancer	PA1338-P4 Pancreatic cancer	PA1360-P6 Pancreatic cancer	
5	PA1383-P7 Pancreatic cancer	PA1443-P7 Pancreatic cancer	PA6236-P11 Pancreatic cancer	PA6285-P3 Pancreatic cancer	PA1290-P2 Pancreatic cancer	PA1322-P4 Pancreatic cancer	PA1390-P3 Pancreatic cancer	PA2021-P2 Pancreatic cancer	PA2022-P5 Pancreatic cancer	PA2031-P8 Pancreatic cancer	PA2425-P5 Pancreatic cancer	PA3546-P3 Pancreatic cancer	
6	PA6284-P4 Pancreatic cancer	PA2847-P12 Pancreatic cancer	PA2924-P4 Pancreatic cancer	PA3023-P6 Pancreatic cancer	PA0372-P5 Pancreatic cancer	PA0787-P5 Pancreatic cancer	PA1168-P4 Pancreatic cancer	PA1170-P2 Pancreatic cancer	PA1194-P2 Pancreatic cancer	PA1198-P2 Pancreatic cancer	PA1221-P3 Pancreatic cancer	PA1222-P3 Pancreatic cancer	
7	PA1229-P2 Pancreatic cancer	PA1252-P2 Pancreatic cancer	PA1280-P3 Pancreatic cancer	PA1301-P4 Pancreatic cancer	PA1332-P2 Pancreatic cancer	PA1337-P3 Pancreatic cancer	PA1338-P4 Pancreatic cancer	PA1360-P6 Pancreatic cancer	PA1383-P7 Pancreatic cancer	PA1443-P7 Pancreatic cancer	PA6236-P11 Pancreatic cancer	PA6285-P3 Pancreatic cancer	
8	PA1290-P2 Pancreatic cancer	PA1322-P4 Pancreatic cancer	PA1390-P3 Pancreatic cancer	PA2021-P2 Pancreatic cancer	PA2022-P5 Pancreatic cancer	PA2031-P8 Pancreatic cancer	PA2425-P5 Pancreatic cancer	PA3546-P3 Pancreatic cancer	PA6284-P4 Pancreatic cancer	PA2847-P12 Pancreatic cancer	PA2924-P4 Pancreatic cancer	PA3023-P6 Pancreatic cancer	
												Mark point	

Map of pancreatic cancer tumor microarray

Please follow standard IHC staining procedure.

Visit <https://www.crownbio.com/databases> to view detailed information for each model. Access to the databases requires online registration via the link provided.