



TMA036

This sarcoma tumor microarray is derived from patient derived xenograft (PDX) tumor samples. Each tumor sample is formalin fixed, paraffin embedded, and spotted in triplicate on the TMA.

Tumor Microarray Profile

Catalog #: TMA-HP-SA-036

Model Type: PDX

Host Species: Mouse

Models On TMA: 30

Cores/Model: 3

Cores On TMA: 90

Core Size: 1 mm diameter

Fixation Method: Formalin fixed, paraffin embedded

Applications: H&E and IHC staining

Storage Conditions: 4°C Stable for 12 months from date of purchase when stored as recommended

Slide Preparation: Slides are sealed with wax; heat slide at 60° for 1 hour before de-waxing

QC: Each TMA slide has >90% core occupancy. H&E and Ki67 IHC staining have been performed on a representative slide of each lot of TMA product.

Tumor Microarray Map

Slice Label	TMA-HP-SA-036											
	A	B	C	D	E	F	G	H	I	J	K	L
1	SA0271-P6 Sarcoma	SA0319-P7 Sarcoma	SA0639-P5 Sarcoma	SA2398-P12 Sarcoma	SA2466-P9 Sarcoma	SA3224-P4 Sarcoma	SA3283-P7 Sarcoma	SA3296-P5 Sarcoma	SA3302-P5 Sarcoma	SA3313-P4 Sarcoma	SA3829-P4 Sarcoma	SA3831-P5 Sarcoma
2	SA3838-P5 Sarcoma	SA3839-P5 Sarcoma	SA3840-P3 Sarcoma	SA3841-P9 Sarcoma	SA3849-P4 Sarcoma	SA4008-P4 Sarcoma	SA4013-P4 Sarcoma	SA4028-P4 Sarcoma	SA4033-P3 Sarcoma	SA4061-P5 Sarcoma	SA4062-P5 Sarcoma	SA4078-P3 Sarcoma
3	SA4094-P4 Sarcoma	SA4098-P3 Sarcoma	SA4109-P2 Sarcoma	SA4117-P2 Sarcoma	SA4134-P5 Sarcoma	SA4136-P4 Sarcoma	SA0271-P6 Sarcoma	SA0319-P7 Sarcoma	SA0639-P5 Sarcoma	SA2398-P12 Sarcoma	SA2466-P9 Sarcoma	SA3224-P4 Sarcoma
4	SA3283-P7 Sarcoma	SA3296-P5 Sarcoma	SA3302-P5 Sarcoma	SA3313-P4 Sarcoma	SA3829-P4 Sarcoma	SA3831-P5 Sarcoma	SA3838-P5 Sarcoma	SA3839-P5 Sarcoma	SA3840-P3 Sarcoma	SA3841-P9 Sarcoma	SA3849-P4 Sarcoma	SA4008-P4 Sarcoma
5	SA4013-P4 Sarcoma	SA4028-P4 Sarcoma	SA4033-P3 Sarcoma	SA4061-P5 Sarcoma	SA4062-P5 Sarcoma	SA4078-P3 Sarcoma	SA4094-P4 Sarcoma	SA4098-P3 Sarcoma	SA4109-P2 Sarcoma	SA4117-P2 Sarcoma	SA4134-P5 Sarcoma	SA4136-P4 Sarcoma
6	SA0271-P6 Sarcoma	SA0319-P7 Sarcoma	SA0639-P5 Sarcoma	SA2398-P12 Sarcoma	SA2466-P9 Sarcoma	SA3224-P4 Sarcoma	SA3283-P7 Sarcoma	SA3296-P5 Sarcoma	SA3302-P5 Sarcoma	SA3313-P4 Sarcoma	SA3829-P4 Sarcoma	SA3831-P5 Sarcoma
7	SA3838-P5 Sarcoma	SA3839-P5 Sarcoma	SA3840-P3 Sarcoma	SA3841-P9 Sarcoma	SA3849-P4 Sarcoma	SA4008-P4 Sarcoma	SA4013-P4 Sarcoma	SA4028-P4 Sarcoma	SA4033-P3 Sarcoma	SA4061-P5 Sarcoma	SA4062-P5 Sarcoma	SA4078-P3 Sarcoma
8	SA4094-P4 Sarcoma	SA4098-P3 Sarcoma	SA4109-P2 Sarcoma	SA4117-P2 Sarcoma	SA4134-P5 Sarcoma	SA4136-P4 Sarcoma						
												Mark point

Map of sarcoma tumor microarray slide

Please follow standard IHC staining procedure.

Visit <https://www.crownbio.com/databases/hubase> to view detailed information for each model. Access to HuBase™ requires online registration via the link provided.